Obelis s.a.  
Bd. General Wahis 53  
1030 Brussels, BELGIUM  
Tel: + (32) 2. 732.59.54  
Fax: + (32) 2. 732.60.03 Email: [mail@obelis.net](mailto:mail@obelis.net)

Distributed by:  
Mizuho America, Inc.  
30057 Ahern Avenue  
Union City, CA 94587  
510-324-4500  
800-699-2547  
[www.Mizuho.com](http://www.Mizuho.com)

Vascular Technology Incorporated

12 Murphy Drive  
Nashua, NH 03062  
**Toll Free** 800.550.0856  
**Outside USA** 603.594.9700  
**Fax** 603.594.0092  
[www.vti-online.com](http://www.vti-online.com)  
**E-mail** [info@vti-online.com](mailto:info@vti-online.com)

## Instructions for Use

### 20 MHz Doppler Systems

REF. 07-150-02, 07-150-20

© 2018, all rights reserved  
201145 Rev 0

## Table of Contents

English.....	ENG-1
Spanish.....	ESP-39
French .....	FRA-77
German.....	DEU-117
Italian .....	ITA-155

# Table of Contents

Technical Description.....	3
Intended Use.....	4
Warnings.....	4
Cautions.....	6
Operation.....	7
Environment.....	7
Power.....	8
Physical – REF 07-150-20.....	8
Physical – REF 07-150-02.....	8
Explanation of Symbols.....	8
Transceiver Description.....	11
System Block Diagram.....	14
Setup.....	15
Battery Installation.....	15
Power Supply Installation.....	15
Transceiver Placement.....	16
Immunity Test.....	17
Performance Criteria.....	26
Essential Performance.....	26
Doppler Probe Selection.....	27
Doppler Probe Connection.....	27
Preparation for Use.....	28
Flow Determination.....	28
Trouble Shooting Guide.....	31
Service.....	34
Maintenance and Cleaning.....	34
Battery Replacement.....	35
Accessories & Parts.....	36
Limited Warranty.....	37

## Technical Description

The 20 MHz Doppler systems are designed, developed and manufactured by Vascular Technology Inc. The system is available as battery – operated or battery and external power supply operated models. The REF 07-150-20 Doppler system is externally DC powered, pulsed Doppler ultrasound systems designed for the evaluation of blood velocity in vessels. The REF 07-150-20 Doppler system can also be powered by eight (8) standard AA (LR6) alkaline batteries. The REF 07-150-02 Doppler system can only be powered by eight (8) standard AA (LR6) alkaline batteries and are also designed for the evaluation through sound of blood velocity in vessels. The system and its probes are and distributed by Mizuho America.

The 20 MHz Doppler systems allow the user to interrogate vessels. A Doppler transducer (probe), which plugs into the Doppler system, emits a pulsed ultrasonic signal. A varying audible signal is produced when the probe is placed upon a vessel within which there is flow. The frequency (i.e., pitch) of the signal is proportional to the blood velocity within the vessel. Distinctive tonal patterns are produced which are indicative of the flow pattern in terms of velocity vs. time. The volume of the tone may be adjusted by means of a control located on the Doppler system. A transmitter in the Doppler system periodically drives the ultrasonic transmitting crystal located at the tip of the probe. The ultrasonic waves generated by the crystal travel through the tissue just under the probe tip in a fairly narrow beam. They are then reflected back towards the probe whenever they encounter a boundary between tissues of different densities. During the intervals when the unit is not transmitting, the probe passes any reflected signals that it receives to a receiving circuit. This circuit amplifies the returning echoes, compares their

frequency to that of the transmitted signal and converts any frequency differences into an audible tone.

## Intended Use

The 20 MHz Doppler systems are intended for the intraoperative and transcutaneous evaluation of blood flow. Indications for Use form filed with the FDA for the Doppler Probes lists the clinical applications as Intraoperative (microvascular and vascular), Intraoperative Neurological, Transesophageal, Transrectal, Laparoscopic and Peripheral Vascular.

### European Union Notice:

The Doppler probes are intended for general use and each probe is intended to be used in multiple surgical specialties. **The Doppler probes are not intended specifically for use in direct contact with the central nervous system (brain, meninges and spinal cord).** The Doppler probes are not intended specifically to control, diagnose, monitor or correct a defect of the heart or of the central circulatory system through direct contact with these parts of the body. The Doppler probes are not intended to be dedicated disposable cardiovascular surgical instruments. The user must follow all Warnings, Cautions and Contraindications associated with this device.

All Doppler probes are designed to be safe and effective when used properly. To reduce the risk of hazards to a level as low as reasonably achievable, and to limit exposure to ultrasound, turn the unit off when not in use.

## Warnings

**WARNING:** *This label alerts the operator to possible bodily injury if procedures are not followed exactly.*

**WARNING:** Never reuse single-use disposable probes. Reuse

may lead to cross contamination and mechanical damage. No proven method exists which can eliminate the possibility of transmitting prion-based brain wasting disease such as variant Creutzfeldt-Jakob Disease (vCJD). Probes which come in contact with brain tissue must be disposed of by incineration.

**WARNING:** Never sterilize the Doppler system with autoclave, ultraviolet, gamma radiation, gas, steam, or heat sterilization techniques. Severe damage and personal injury could result.

**WARNING:** There are no user serviceable components inside this device other than the batteries. Disassembly of the internal components of this unit may result in circuit damage. All servicing should be referred to the factory or Mizuho America. Do not modify this device.

**WARNING:** Not for use in OXYGEN ENRICHED atmospheres.

**WARNING:** Remove batteries during prolonged storage. Do not install backwards, charge, put in fire, or mix with other battery types. Batteries may explode or leak causing injury. Replace all batteries at the same time.

**WARNING:** Equipment not suitable for use in the presence of a FLAMMABLE ANESTHETIC MIXTURE WITH AIR OR WITH OXYGEN OR NITROUS OXIDE.

**WARNING:** The Doppler system should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Doppler system should be observed to verify normal operation in the configuration in which it will be used.

## Cautions:

**CAUTION:** *This label alerts the operator to possible equipment or software damage if procedures are not followed correctly.*

**CAUTION:** This Doppler device is not intended for fetal use.

**CAUTION:** Prior to use, inspect probe for damage, such as cracks that could allow for the ingress of fluids into the probe, and sharp edges.

**CAUTION:** Use only AA (LR6) alkaline batteries.

**CAUTION:** Use Doppler system only with compatible power supply.

**CAUTION:** Because the Doppler system needs to be sensitive to very weak signals from blood flow, by design it may be susceptible to picking up interference through the coaxial cable that connects the Doppler probe to the Doppler system.

**CAUTION:** Use Doppler system only with compatible Doppler probes.

**CAUTION:** The Doppler probe is delicate. Do not drop or strike against hard surfaces. Avoid excessive mechanical pressure on the probe or excessive tension on the probe cable as the cable may break or cause the sensor to pull back from and not function. Check to ensure connectors are not loose.

**CAUTION:** To avoid biological hazards, properly dispose of probe according to local regulations.

**CAUTION:** Properly dispose of batteries according to local regulations. The Doppler system and power supply may be returned to manufacturer for proper disposal.

**CAUTION:** The Doppler system should not be used in the presence of any high frequency equipment, including high frequency surgical generators.

## Operation

Transmission Frequency	20 MHz
Transmission Characteristic	Pulsed transmission, continuous reception


## Environment

Ambient Operating Temperature Range	+15° to +40°C.
Ambient Operating Humidity Range	30% to +75% RH
Ambient Shipping Humidity Range	30% to +75% RH
Ambient Shipping Temperature Range	-15° to +70° C.
Ambient Storage Temperature Range	-15° to +40° C.
Ambient Storage Humidity Range	30% to +75% RH
IPX 0 (Transceiver)	No special protection
IPX 7 (Probe)	Protected against water immersion - Immersion for 30 minutes at a depth of 1 meter.
Surface Temperature – probe	Less than 43°C

Power	
REF 07-150-20 Mizuho Doppler System	8 AA (LR6) alkaline batteries and External Power Source, A/C to D/C Power Supply
REF 07-150-02 Mizuho Doppler System	8 AA (LR6) alkaline batteries
Power Requirements	12 VDC

Physical – REF 07-150-20	
Dimensions	6.5 in. D X 10 in. W X 4 in. H. (165 mm X 254 mm X 101 mm), nom
Weight	2.6 lb, (1.18 kg), nom.
Physical – REF 07-150-02	
Dimensions	4.25" D X 3.75" W X 5.25" H (105mm X 95mm X 135mm), nom
Weight	1.4 lb, (0.6 kg) nom

Explanation of Symbols	
Sterilized by Ethylene Oxide	

Do Not Re-use	
Follow Instructions for Use	
Date Manufactured	
Use By	
Type CF Applied Part	
Catalog Number	
Lot	
Serial Number	
Keep Dry	
Direct Current	
Protection From Ingress of Water	
RF transmitter	
Authorized representative in the European Community	
Limit of temperatures	

Do Not Re-sterilize

# Doppler System Description

## REF 07-150-20


**Fig. 4-1** REF 07-150-20 Doppler Transceiver (20 MHz System is depicted)

1. Power Switch: A push-button control that when depressed turns the unit ON. Power ON is indicated by illumination of one of the four green Volume Setting LEDs. When the push-button is depressed a second time the unit is turned OFF.

2. Volume Indicator: A series of green LEDs that indicates the volume of the audible Doppler signal.
3. Volume Increase Switch: A push-button switch that when depressed and held will increase the volume of the audible Doppler signal.
4. Volume Decrease Switch: A push-button switch that when depressed and held will decrease the volume of the audible Doppler signal.
5. Low Battery Indicator: A yellow LED which, when illuminated, indicates that the batteries are nearing the end of their useful life. The system will automatically shut OFF when the battery voltage is too low to maintain proper operation of the unit.

# Doppler System Description

## REF 07-150-02

**Fig. 4-2 REF 07-150-02 Doppler System**

1. Power Switch: A push-button control that when depressed turns the unit ON. Power ON is indicated by illumination of any of the four green Volume Setting LED's. When the push-button is depressed a second time the unit is turned OFF.
2. Battery Low Indicator: A yellow LED which, when illuminated, indicates that the batteries are nearing the end of their useful life. The system will automatically shut OFF when the battery voltage is too low to maintain proper operation of the unit.
3. Volume Indicator: A series of green LED's that indicates the


volume of the audible Doppler signal.

4. Volume Increase Switch: A push-button switch that when depressed and held will increase the volume of the audible Doppler signal.

5. Volume Decrease Switch: A push-button switch that when depressed and held will decrease the volume of the audible Doppler signal.

### System Block Diagram

**Fig. 4-2 System Block Diagram**

**CAUTION:** Prior to use, inspect probe for damage, such as cracks that could allow for the ingress of fluids into the probe, and sharp edges.

Carefully unpack your Doppler system. Inspect the Doppler system for damage. If the Doppler system is missing or any damage is found, contact Mizuho America for further instructions.

## Setup

### Battery Installation

**CAUTION:** Use only AA (LR6) alkaline batteries.

Install the batteries per the instructions in the Service Information section of this manual.

### Power Supply Installation

**CAUTION:** On models that use a power supply, use only with compatible power supply.

On models that use a power supply, connect power supply to wall grounded hospital grade outlet using hospital grade cord set. Connect power connector to D/C outlet on unit.

On models that use a power supply, patient isolation from the mains is accomplished in the following ways: First, an external medical grade power supply is used to provide D.C. operating voltage for the Doppler unit. Both D.C. power output lines from the power supply are isolated from the mains in the D.C. power

supply. There is no connection between the "green" safety ground and the Doppler unit. Additionally, Doppler probes that make patient contact are further isolated by use of an RF isolation transformer that isolates both signal and ground leads between the Doppler unit and the Doppler probe. The transducer coaxial shield is connected to the Doppler unit chassis ground through two safety capacitors. Both the capacitors and the isolation transformer have a 4000VAC dielectric withstand voltage. The final isolation mechanism is the cable insulation and potting of the probe providing an additional insulation layer between the isolated electrical signals and the patient.

## Doppler System Placement

**CAUTION:** Because the Doppler system needs to be sensitive to very weak signals from blood flow, by design, it may be susceptible to picking up interference through the coaxial cable that connects the Doppler probe to the Doppler system.

The Doppler system needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information located in this Instructions for Use manual.

Portable and mobile RF communications equipment can affect the Doppler system.

**WARNING:** The Doppler system should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Doppler system should be observed to verify normal operation in the configuration in which it will be used.

Place the Doppler System on a suitable stand, cart or table outside the sterile field, just behind the physician who will be using the Doppler system. The Doppler system should not be used in the presence of any high frequency equipment, including high frequency surgical generators. The transceiver has a rating of IPX-0. Keep the transceiver away from all open liquids.

### Immunity Test

<b>Table 1 Guidance and manufacturer’s declaration – electromagnetic emissions</b>		
The Doppler systems are intended for use in the electromagnetic environment specified below. The user of the Doppler system should assure that it is used in such an environment.		
<b>Emissions Test</b>	<b>Compliance</b>	<b>Electromagnetic Environment -- guidance</b>
RF Emissions, CISPR 11	Group 1	The Doppler system uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	The Doppler systems meet the conducted and radiated performance requirements for non-life supporting equipment and also meet the harmonic emissions, voltage dips and variations and
Power Harmonic emissions	Not Applicable (see Note)	

<p>Voltage fluctuations/flicker emissions</p>	<p>Not Applicable (see Note)</p>	<p>voltage fluctuation (flicker) requirements for non-life supporting equipment pursuant to IEC 60601-1-2:2007 and CISPR 11, A1 &amp; A2, and IEC 61000-3-3.</p> <p>The Doppler systems are suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded:</p> <p>Warning: The Doppler systems are intended for use by healthcare professionals only. The Doppler systems may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the Doppler system or shielding the location.</p>
<p>NOTE: Outside the scope of the standard.</p>		

**Table 2 Guidance and manufacturer's declaration –  
electromagnetic Immunity**

The Doppler systems are intended for use in the electromagnetic environment specified below. The customer or the user of the Doppler System should assure that it is used in such an environment.

<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment -- guidance</b>
Electrostatic discharge (ESD)  IEC 61000-4-2	+/- 6 kV contact  +/- 8 kV air	+/- 2, +/- 4 and +/- 6 kV Contact Discharge.  +/- 2, +/- 4, and +/- 8 kV Air Discharge	Floors should be wood, concrete or ceramic tile.  If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst  IEC 61000-4-4	+/- 2 kV AC Mains  +/- 1 kV I/O Ports	+/- 2 kV AC Mains  +/- 2 kV I/O Ports	Mains power quality should be that of a typical commercial or hospital environment.

<p>Surge</p> <p>IEC 61000-4-5</p>	<p>+/- 1 kV line(s) to line(s)</p> <p>+/- 2 kV line(s) to earth</p>	<p>+/- 1 kV line to line (DM)</p> <p>+/- 2 kV line(s) to earth (CM)</p>	<p>Mains power quality should be that of a typical commercial or hospital environment.</p>
<p>Voltage dips, short interruptions and voltage variations on power supply input lines.</p> <p>IEC 61000-4-11</p>	<p>&lt;5% <math>U_T</math> (&gt;95 % dip in <math>U_T</math>) for 0,5 cycle.)</p> <p>40 % <math>U_T</math> (60 % dip in <math>U_T</math>) for 5 cycles</p> <p>70 % <math>U_T</math> (30 % dip in <math>U_T</math>) for 25 cycles</p> <p>&lt;5 %</p>	<p>10 ms &gt;95 % dip (0.5 period)</p> <p>100 ms 60 % dip (5 periods)</p> <p>500 ms 30 % dip (25 periods)</p> <p>5 seconds &gt; 95 % dip/interruptions (300 periods)</p>	<p>Mains power quality should be that of a typical commercial or hospital environment.</p> <p>If the user of the Doppler System requires continued operation during power mains interruptions, it is recommended that the Doppler System be powered from an uninterruptible power supply or the built-in battery.</p>

	U <sub>T</sub> (>95 % dip in U <sub>T</sub> ) for 5 s		
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m 50 Hz & 60 Hz three orthogonal orientations	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U <sub>T</sub> is the A.C. mains voltage prior to application of the test level.			

**Table 3 Guidance and manufacturer's declaration –  
electromagnetic immunity**

The Doppler systems are intended for use in the electromagnetic environment specified below. The customer or the user of the Doppler System should ensure that it is used in such an environment.

IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT -- GUIDANCE
<p>Conducted RF IEC 61000-4-6</p>	<p>3 Vrms 150 kHz to 80 MHz □</p>	<p><math>[V_1]V = 3</math> Vrms □</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Doppler system, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance:</b></p> $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P} = 1.2\sqrt{P}$

<p>Radiated RF IEC 61000- 4-3</p>	<p>3 V/m 80 MHz to 2,5 GHz</p> <p style="text-align: center;">□</p>	<p><math>E_1 = 3 \text{ V/m}</math></p> <p style="text-align: center;">□</p> <p style="text-align: center;">□</p>	$d = \left[ \frac{3.5}{E_1} \right] \sqrt{P} = 1.2 \sqrt{P}$ <p>80 MHz to 800 MHz</p> $d = \left[ \frac{7}{E_1} \right] \sqrt{P} = 2.3 \sqrt{P}$ <p>800 MHz to 2.5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m)</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: ((⚡))</p>
---	---	---	--

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Doppler system is used exceeds the applicable RF compliance level above, the Doppler system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Doppler system.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Table 4 Recommended separation distances between portable/mobile RF communications equipment and the Doppler Systems.**

The Doppler systems are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Doppler system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Doppler system as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter  w	Separation distance according to frequency and power of transmitter  (m)		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2,5 GHz
	$d = 1.2\sqrt{P}$	$d = 1.2\sqrt{P}$	$d = 2.3\sqrt{P}$
0,01 □	0,12 □	0,12 □	0,23
0,1	0,38	0,38	0,74
1	1,2	1,2	2,3
10	3,8	3,8	7,4
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter power in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## Performance Criteria

Failures include any time the unit does not produce an audible signal when detectable flow is present. In addition to component malfunction, failures also include units that produce a false audible signal that is indistinguishable from a signal produced by flow. Non-intentional audible signal tones are allowed to be produced by the unit, so long as they cannot be easily mistaken for flow.

## Essential Performance

The Doppler system is not intended to be used as the sole indicator of blood flow detection. Other clinical means should be utilized to check for failure or degradation of the Doppler system.

The equipment or system may exhibit degradation of performance (e.g. Deviation from specifications) that does not affect essential performance or safety. If the Doppler probe cable is in close proximity to electrically active conductors, such as an electro-surgery cables or an electronic equipment chassis, signals from the cables or chassis may be picked up by the Doppler system and produce audible signals that will

be heard along with the desired Doppler audio. This interference is easily distinguished from blood flow and is remedied by moving the transducer cable away from the source of the interference. Always recheck the Doppler System for proper function after this interference has occurred.

**WARNING:** Equipment not suitable for use in the presence of a FLAMMABLE ANESTHETIC MIXTURE WITH AIR OR WITH OXYGEN OR NITROUS OXIDE.

**WARNING:** Not for use in OXYGEN ENRICHED atmospheres.

## Doppler Probe Selection

**CAUTION:** Use Doppler system only with compatible Doppler probes.

Doppler systems are designed to function only with Doppler probes that are compatible. Patient contact areas of the probe are rated at IPX-7. Keep connectors away from all liquids. Acoustical output tables and information required by IEC60601-2-37 can be found in the inserts provided with the Doppler probes.

## Doppler Probe Connection

Using sterile technique when appropriate, remove the sterile Doppler probe from its packaging. Hand-off the probe's connector to someone outside the sterile field. Attach the probe's connector to the coaxial receptacle on the transceiver front panel.

**WARNING:** Never sterilize the Doppler system with autoclave, ultraviolet, gamma radiation, gas, steam, or heat sterilization techniques. Severe damage and personal injury could result.

## Preparation for Use

Plug the power cord into a hospital grade electrical outlet or allow the unit to be powered from the internal battery source. Turn the transceiver on by depressing the Power Switch. The initial volume level will be maintained from the previous use.

Adjust the volume by depressing and holding the Volume Increase or Decrease push-button to an intermediate level. Some “white” noise (white noise sounds similar to a radio that is tuned between stations) may be heard from the Doppler system’s speaker. To verify that the system is operational, gently draw the tip of the Doppler probe, using sterile technique, along any convenient sterile surface. This will produce a fairly loud rasping noise, confirming that the system is operational.

**CAUTION:** Prior to use, inspect probe for damage, such as cracks that could allow for the ingress of fluids into the probe, and sharp edges

**CAUTION:** The Doppler probe is delicate. Do not drop or strike against hard surfaces. Avoid excessive mechanical pressure on the probe or excessive tension on the probe cable. Check to insure connectors are not loose.

## Flow Determination

Place the tip of the probe directly on the vessel or other site to be evaluated, orienting the probe as shown in figure 6-1. Adjust the

angle between the probe and the vessel until the maximum audible signal is obtained. Adjust the volume control on the Doppler system to the desired level. If any flow is detected, the pitch of the resultant audible signal will correspond to its velocity, with higher pitches indicating higher velocities. The probe may be moved to various sites as required.

When the entire procedure is finished, turn the Doppler system's power OFF by depressing the Power switch. Clean the Doppler system, if necessary, as described in the Service Information section. Properly dispose of probe according to local regulations.

**CAUTION:** To avoid biological hazards, properly dispose of probe according to local regulations.

**WARNING:** Never reuse single-use disposable probes. Reuse may lead to cross contamination and mechanical damage. No proven method exists which can eliminate the possibility of transmitting prion-based brain wasting disease such as variant Creutzfeldt-Jakob Disease (vCJD). Probes which come in contact with brain tissue must be disposed of by incineration.


Fig. 6-1 Representative Doppler Probe Orientation – Refer to Instructions for Use supplied with individual probes for specific details.

## Trouble Shooting Guide

Symptoms	Possible Problems & Solutions
Weak sound output, even at maximum volume setting.	The flow that is being heard is somewhat deeper than this unit is designed to detect.
	Weak batteries (Yellow LED). Replace.
	Doppler probe may be defective. Return probe; contact factory.
	Doppler system may be defective. Contact factory.

<p>“White” noise occurs at maximum volume setting and drawing the probe tip over a surface results in a rasping noise, but probe does not detect flow.</p>	<p>Doppler probe not connected. Connect.</p> <hr/> <p>Doppler probe is correctly evaluating a zero velocity condition. No problem.</p> <hr/> <p>Doppler probe is not positioned correctly. Review Operation section.</p> <hr/> <p>Doppler probe may be defective. Return to the factory for repair.</p> <hr/> <p>Doppler system may be defective. Return for repair.</p>
--	--

<p>“White” noise occurs at maximum volume setting, but drawing the probe tip over a surface does not result in a rasping noise.</p>	<p>Doppler probe not connected. Connect.</p> <p>.....</p> <p>Doppler probe may be defective. Contact factory.</p> <p>.....</p> <p>Doppler system may be defective. Contact factory.</p>
<p>No sound whatsoever, at any volume control setting; low battery indicator not illuminated.</p>	<p>Power Supply not connected. Connect.</p> <p>.....</p> <p>Power Supply may be defective. Contact factory.</p> <p>.....</p> <p>Transceiver may be defective. Contact factory.</p>

If the problem cannot be corrected after making the above checks and adjustments, call Mizuho America for additional help or return authorization at (800-699-CLIP (2547) between the hours of 9:00 A.M. and 5:00 P.M. Pacific Standard Time, Monday through Friday.

Pacific Time, Monday through Friday.

## Service

If you have any problems operating your instrument, read the instructions again, paying particular attention to the sequential order of steps; or consult the Troubleshooting section.

**WARNING:** There are no user serviceable components inside this device other than the batteries. Disassembly of the internal components of this unit may result in circuit damage. All servicing should be referred to the factory.

## Maintenance and Cleaning

The Doppler system requires little maintenance. Keep it clean and free of dust. The exterior may be cleaned using the following steps:

1. After every use, check the Doppler system for any sign of damage or wear.
2. Wipe the Doppler system with a dry or water-moistened soft cloth.
3. Wipe the Doppler system with Isopropyl alcohol (70%) moistened soft cloth or wipe. Do not pour 70% isopropyl alcohol directly on the Transceiver. Always use a moistened soft cloth or wipe. Allow to air dry before use.
4. Check the Doppler system for any residual organic material. If any is present, remove it and disinfect the Doppler system again.

The Doppler system should not contact mucus membranes, blood, or compromised tissue, and is not used in sterile fields.

**WARNING:** Never sterilize the Doppler system with autoclave, ultraviolet, gamma radiation, gas, steam, or heat sterilization techniques. Severe damage and personal injury could result.

## Battery Replacement

Make sure that the Doppler system's power is turned OFF by depressing the Power Switch push-button so that no green LED's are illuminated.

Remove the battery cover door of the Doppler system, exposing the battery compartment. Slide the battery holder out of the compartment. Remove the old batteries and replace with fresh, AA (LR6) size alkaline batteries only, taking care to observe the battery polarity markings that are molded into the battery holder. Replace the battery holder and battery door cover of the transceiver. Properly dispose of batteries according to local regulations.

**WARNING:** Remove batteries during prolonged storage. Do not install backwards, charge, put in fire, or mix with other battery types. May explode or leak causing injury. Replace all batteries at the same time.

## Accessories & Parts

The use of ACCESSORIES, transducers and cables other than those specified, with the exception of transducers and cables sold by the MANUFACTURER of the ME EQUIPMENT or ME SYSTEM as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of the ME EQUIPMENT or ME SYSTEM.

Item	Catalog Number
Power Supply	108110-SUPPLY
Power Cord	108110-US
Batteries	101325
Instructions for Use	201083

<b>Probes for use with 07-150-02 and 07-150-20</b>		
Item	Catalog Number	Type
20 MHz Curved, Single Use Transducer	07-150-22	CF
20 MHz Single Use Probe	07-150-20	CF
20 MHz Single Use Probe	108380	CF
20 MHz Bayonet, Single Use Probe	07-150-07	CF
20 MHz Bayonet, Slim, Single Use	07-150-10	CF

Probe		
20 MHz Bayonet, Mini, Slim, Single Use Probe	07-150-12	CF

## Limited Warranty

The Doppler system is warranted for one (1) year from the date of shipment from the factory against defects in materials and workmanship. Defective Doppler systems will be repaired or replaced, at VTI's option, when returned prepaid to the factory within this year.

The customer assumes full responsibility that this equipment meets the specifications, capabilities and other requirements of the customer. VTI makes no warranty of fitness for a particular purpose except as provided herein.

The customer assumes full responsibility for the proper installation, operation and maintenance of this equipment as described in this manual as well as other instructions that may be provided by VTI. This warranty is void if the equipment has been mishandled, operated outside of its specified operating or environmental limits or otherwise subjected to improper or abnormal use.

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