MIZUHO

CE
0123

SUGITA

Titanium ANEURYSM CLIPS II

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Emergo Europe B.V.
Prinsessegracht 20, 2514 AP, The Hague
The Netherlands

MIZUHO Corporation

3-30-13 Hongo, Bunkyo-ku, Tokyo 113-0033, Japan

TEL(03)3815-3096 FAX(03)3813-5068

<http://www.mizuho.co.jp>

CK17-001-01WW

2019-05 Ver.18

INDEX

English	1	Italiano	67
Deutsche	23	Español	91
Français	45	日本語	113

INSTRUCTIONS FOR THE USE OF SUGITA TITANIUM ANEURYSM CLIPS II

Index

- | | |
|---|---------------------------------|
| 1 . Specifications | 11. Use of MRI |
| 2 . Types | 12. Caution |
| 3 . Spring Type | 13. Defects and Adverse Events |
| 4 . Action Mechanism | 14. Storage - Care and Handling |
| 5 . Surface Finish | 15. Validated Cleaning and |
| 6 . Blade Cross-sectional Shape | Sterilization Procedure |
| 7 . Indication for Use | |
| 8 . Directions for Use | |
| 9 . Contraindications and
Prohibitions | |
| 10. Warning | |

1. Specifications

Sugita Titanium Aneurysm Clips II series offer the permanent type and temporary type.

“Sugita Titanium Aneurysm Clips II” herein referred to as “clip”.

This device is made of titanium 6-aluminium 4-vanadium alloy (Ti-6Al-4V ELI; as specified by ISO 5832-3/ASTM F136) and unalloyed titanium (specified by ISO 5832-2/ASTM F67 G2).

This clip designed for **single use only** is exclusively indicated only for use with the Sugita Titanium Aneurysm Clips II Applying Forceps (hereinafter this is called “clip applier”), and requires cleaning and sterilization prior to use. Reuse may cause the breakage of clips and lead to infection.

ASTM: American Society for Testing and Materials


ISO: International Organization for Standardization

The original design adopts wide opening of the blades providing optimum visibility of the cerebral aneurysms

during surgery.

In order to ensure the traceability, a serial number is permanently marked on each clip.

The closing force of each clip is measured during the final quality assurance stage with a computerized electronic scale at a measuring point located $1/3$ L from the distal end of the blade length (L) of the clip [**Figure 1**] while its blades are opened at a 1.0 mm width. The result of this measurement is recorded in newtons and grams and printed on the device's packaging label.

[Figure 1]

2. Types

Permanent Type

Type	Color
Standard Type	Spring-green
	Blade-silver
Mini Type	Spring-pink
	Blade-silver
Long Type	Spring-blue
	Blade-silver

Temporary Type

Type	Color
Standard Type	Spring-green
	Blade-yellow
Mini Type	Spring-pink
	Blade-yellow

For specific dimensions and shape options, please see this device's catalog.

3. Spring Type : Integral Coil

4. Action Mechanism : Alligator Action

5. Surface Finish

- Clip Body : Satin

- Material Treatment : Anodizing Treatment

6. Blade Cross-sectional Shape: Rolled

7. Indication for Use

The permanent type is indicated for the permanent occlusion of the neck of cerebral aneurysms.

The temporary type is indicated for the temporary occlusion of the neck of cerebral aneurysms.

For the US market

The permanent type is indicated for the permanent

occlusion of the neck of cerebral aneurysms and peripheral blood vessels.

The temporary type is indicated for the temporary occlusion of the neck of cerebral aneurysms and peripheral blood vessels.

8. Directions for Use


This device is packaged non-sterile and it must be properly cleaned and sterilized before use.

Each clip must be held with the specified clip applicator.

Each clip applicator is etched to read; e.g. “FOR (type of clip) Only”.

When holding the clip, it must be verified that the spring of clip is correctly placed in the holding jaws of the clip applicator as shown in [Figure 2].

Improper holding of the clip may cause damage or ejection from the clip applicator. As this may present a risk during surgery, correct holding of the clip must always be kept in mind.

[Figure 2]

- Verify correct placement of the clip in the recesses of the holding jaws of clip applicator.
- Do not hold the clip in angled positions.
- If the spring of the clip is not properly held with the holding jaws of clip applicator, the clip may slip or flick off.

This may cause damage to the spring section of the clip with possible fracture of the spring.

It is very important to verify correct placement of this device to the target site.

It is recommended that the hospital record the type, lot and serial numbers of the clip used for each patient. Patient/Trace labels are included in the clip package for this purpose.

9. Contraindications and Prohibitions

1) Prohibition of reuse

This device is for single-use only.

Once the blades of the clip have been opened or the clip has come into contact with blood or tissue of a patient, it must be discarded.

Reuse may cause the breakage of clips and lead to infection. [Cause of breakage/risk of infection]

2) Prohibition of use for unintended purposes

This device should not be used for purposes other than indication for use.

3) Use only the exclusive clip applier when handling this clip to avoid the changes in the closing forces or mechanical features of the clip. Do not use other

manufacturers' clip applicators with this clip as there is risk of harm to the patient, breakage or damage to the clip.

4) Do not open the blades of clip directly with fingers. This may cause corrosion or alter the proper holding force. The blades of clip should be opened only when the spring of clip is placed in a proper position of the clip applicator.

5) This device should not be bent or remodeled.

6) In order to preclude any possible Galvanic reaction, do not use this device with implants of other manufacturers or with implants made of different materials. These Galvanic reactions may cause corrosion or alter the proper holding force.

7) This device must not be used in patients with serious allergies to metals or exogenous materials.

8) This device must not be used in patients with infectious diseases or sepsis as a complication.

9) This device must not be sterilized using low-temperature hydrogen peroxide gas plasma.

This product is not compatible with sterilization using low-temperature hydrogen peroxide gas plasma.

[Sterilization using such method may cause surface color deterioration and affect its characteristics.]

10. Warning

1) Full explanation should be given to each patient about the use of this device as well as the benefits/risks of surgery involving its use.

2) Closing status of the blade of clip should be verified immediately following surgery.

[The most serious adverse events expected to be reported in conjunction with the aneurysm clip may include displacement of the clip, slip-out of the clip, ejection of the clip from clip applicator, clip fracture, laceration of the aneurysm, cerebrovascular injury, and hemorrhage. In addition, in consideration of the structure of the clip, risk of failures/adverse events may be increased by applying only the tips of clips, especially the ones with long

blades, to the target tissue.]

3) Prohibition of secondary processing on this product

Secondary processing on this product in the form of, for example, incision by applying impact or vibration on it should not be done. [It may cause breakage.]

4) Prohibition of rough handling

Rough handling of this product such as handling that may cause deformation or damage should be avoided. [It will markedly shorten the life of the implant.]

11. Use of MRI

MR Conditional

Non-clinical testing demonstrated that this device is MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5 tesla or 3 tesla, only
- Maximum spatial gradient magnetic field of 4,000 gauss/cm (40-T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode

Under the scan conditions defined, this device is expected to produce a maximum temperature rise of 2.2°C after 15-minutes of continuous scanning (i.e., per

pulse sequence).

In non-clinical testing, the image artifact caused by this device extends approximately 5-mm from this device when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

NOTE: Do not take the clip applier into the MR environment. This clip applier contains material which has been determined to be MR Unsafe.

12. Caution

Use of this device is restricted to surgeons with proper training and experience in the use of intracranial aneurysm clips.

13. Defects and Adverse Events

1) Serious adverse events

- Displacement of the clip (slip-out)
- Ejection of the clip from clip applicator
- Breakage of the clip
- Rupture of a cerebral aneurysm
- Cerebrovascular injury
- Hemorrhage
- Death

2) Other adverse events

- Infection of surgical wound
- General surgical complications
- Delayed clip displacement due to the application of magnetic field
- Occlusion of adjacent blood vessels
- Foreign body reaction

14. Storage - Care and Handling

This device is prepared in a non-sterile package and should be stored in its package in a cool, dry, and secure location.

Do not use any clip applicator other than the exclusive clip applicator.

15. Validated Cleaning and Sterilization Procedure

Proper cleaning, sterilization and handling of this device is important for assuring that the clips meet sterile specifications prior to surgery as an aid in minimizing the occurrence of infection in the surgical patient.

Other procedures are also important in minimizing contamination in an operating room including the use of a surgical unit with laminar flow equipment at the surgeon's discretion.

This device is supplied nonsterile.

Therefore, this device must be sterilized by the users (healthcare professionals) prior to surgery using validated sterilization processes in accordance with guidelines stipulated by the government.

This clip and instruments must be cleaned to reduce bioburden before first sterilization and re-sterilization. Do not use chlorine-based solutions when cleaning this product. This could cause metallic corrosion.

1) Cleaning

Put the clip in the designated Sugita Clip Case or Sugita Clip Tray (hereinafter this is called "clip case or tray") and close the lid.

<In case of using Ultrasonic cleaner>

(1) Fully immerse the case (tray) in the water (tap water) of the washing tank.

(2) Ultrasonic cleaning

Cleaning time: 15 minutes or more

Detergent: Neutral enzyme type detergent (Diluted concentration 1.0%)

(3) Rinsing

- Make sump water with running water, immerse the case (tray) completely in water (tap water or deionized water), shake the case (tray) 5 times up and down.
- Renew the water and shake the case (tray) 5 times up and down in the water again.
- Perform a series of rinsing operations three or more times in total.

(4) Dry thoroughly.

Recommended drying condition: 110 °C, 30 minutes or more

<In case of using washer disinfector (automatic washing)>

Stage	Procedure	Temperature [°C]	Time [min]	Quality of water	Remarks
I	Preclean	<25	3	Tap water	–
II	Washing	50	10	Tap water	Alkaline detergent (such as cleaning agent for aluminum) that does not affect anodic oxidation coatings Diluted concentration 0.5%
III	Rinsing	<25	4	Deionized water	Set rinsing process program as "1 time (1 minute) × 4 times = 4 minutes"
IV	Hot water disinfection	–	–	Deionized water	Follow the program of the device.
V	Drying	–	–	–	

2) Inspection

Check products visually one by one and do not use if the following abnormalities are found.

- Damage
- Deformation
- Corruption
- Inappropriate blade position (deviation, clearance, tip unevenness, etc.)
- Remaining dirt

Since the color of the surface may be faded in repetitive processing, do not use products that cannot be clearly distinguished with color coding.

3) Packaging

Place clip in the designated clip case or tray and close the lid.

Then pack the case (tray) containing the product using standard medical packaging material (non-woven fabric or sterilized bag) with protective sterile barrier complying with ISO 11607 for steam sterilization.

4) Sterilization

Perform sterilization treatment under sterilization condition "134°C, 3 minutes or more" with a prevacuum type autoclave sterilizer which was confirmed to function properly.

Do not use low temperature hydrogen peroxide gas plasma sterilization. The surface fades and affects its properties.

5) Storage

Do not store the device in high temperatures, or in areas with high humidity and/or where the temperature varies. Please store the device in clean dust-free and dry area.

- General Requirements -

- Handling of this device should be minimized to reduce contamination or damage.
- Sterilization equipment should be validated to assure that it sterilizes specific load configurations properly. The end user is responsible to validate the sterilization process in the individual sterilizer and assure the sterility.
- Sterilization equipment should be cleaned and maintained in accordance with standardized operating procedures specified by its manufacturer.

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