EN

## Instructions for Use - Suction Tube Instruments

### Intended Use

This device is intended for use by a surgeon with professional knowledge of human physiology, in a clinical environment for the purpose of removing fluid and debris during surgery.

### CAUTION

**Federal law restricts this device to sale by or on order of a physician.**

This IFU applies to multiple types of Suction Tube Instruments. Prior to use of this product ensure the instrument is only used by surgeons with the required training, knowledge and/or experience. No additional training is required as Surgeons have an intimate knowledge of anatomy and physiology, and are trained in the use of similar devices as part of their surgical practice. Read, follow and retain these instructions and only use this device for its intended purpose.

### Instruction for Use

Prior to and after each use, inspect to ensure proper function, inspect for damage, and clean and sterilize in accordance with these instructions. If any problems are found related to the device, do not use the device and contact your Mizuho America Representative for support.

### Warnings

- There is a risk of injury and/or failure if the product is not inspected prior to and after each use.
- There is a risk of injury when using this device outside of the field of view.
- There is risk of injury when using this device due to high suction. Prior to use, the pressure flowrate should be set by the user according to the hospital and/or manufacturer's requirements.
- Product should be used according to these instructions and with close attention to patient safety. Failure to follow these and/or the healthcare facility's requirements may cause harm.
- Do not reuse the device when it is used on a patient with, or suspected of having, Creutzfeldt-Jakob Disease (CJD) or variant CJD (vCJD).
- The device must be inspected, cleaned, and sterilized before each surgical procedure.

### Contraindication/Precautions

- The improper use of this device during the handling, surgical use or reprocessing procedures for which it is indicated, may result in damage or device failure.
- To avoid scratching and discoloration, keep ebonized devices separate from other devices.
- Do not use the device in such a manner as to cause excessive stress or strain on the device as misuse may result in misalignment or failure of the device and/or its connections.
- Washing with an improper detergent may result in discoloration or corrosion.
- This device is provided non-sterile.
- Material: Stainless Steel.

### Maintenance, Inspection and Functional Testing

- Prior to and after each use, and after each cleaning, visually inspect all surfaces for damage, distortion, cracks or corrosion and ensure the instrument is functional.
- Visually inspect for damage and wear: device must be free of nicks and

burs; devices with threads must align correctly and tubing should not be kinked or occluded.

- Do not use devices with joints without applying a surgical grade instrument milk/anticorrosive lubricant onto its mating parts as galling could occur and may cause damage to the instrument.
- Do not use if device is defective or damaged.

### Cleaning & Sterilization

#### Decontamination, Cleaning & Drying

- Reprocess all devices as soon as they are ready for cleaning; delaying the time for cleaning may damage or stain the device.
- Dried or affixed surgical residues may make cleaning more difficult or ineffective and may lead to corrosion.
- Residues containing chlorine or chlorides (e.g. from surgical residues, drugs, saline solutions, and water for cleaning, disinfection, and sterilization) may cause corrosion and/or surface damage to stainless steel making the product(s) unusable or unrepairable.
- Remove gross contaminants with a steady stream of lukewarm/cool water, not to exceed 35°C (95°F).
- Do not soak devices in hot water, alcohol, disinfectants or antiseptic solutions. Do not exceed 2 hours soaking in any solution.
- Do not use steel wool, wire brushes or abrasive detergents.
- Devices manufactured from different metals should be processed separately to avoid electrolytic action between the metals resulting in pitting and rusting of the devices.
- Disassemble devices for cleaning. Give special attention to connections and the inside and outside of the suction tubing.
- Where available, use appropriate cleaning accessories to flush suction tubing.
- The product must be rinsed thoroughly with demineralized water prior to washing.
- Follow instructions and warnings as issued by manufacturers of any decontaminants, disinfectants, and cleaning agents used. Wherever possible avoid use of mineral acids or harsh, abrasive agents.
- Use only medical detergents to clean this device.
- Do not use high acid (pH 4.0 or lower) or high alkaline (pH 10 or higher) products for disinfection and cleaning. Neutral pH detergents of 7.0-9.0pH are recommended.
- After washing, rinse each device thoroughly. Do not use saline or chlorinated solutions. Only use distilled or deionized water.

#### Cleaning & Drying: Automated

- Use only validated washer-disinfector machines and low-foaming, non-ionizing cleaning agents and detergents for surgical instruments used in a healthcare environment.
- Follow the manufacturer's instructions for use, which includes warnings, concentrations and recommended cycles, and carefully load the devices.
- Place heavy instruments on the bottom containers, taking care not to place on delicate devices or overload wash baskets.
- Place instruments with curved surfaces facing down to prevent pooling of water.
- Ensure that soft, high purified water that is controlled for bacterial endotoxins is used in the final rinse stage.
- Devices must be dry and all residual moisture removed before sterilization. Use a soft absorbent towel or cloth to dry external surfaces.
- After washing, apply a surgical grade instrument milk/anticorrosive lubricant prior to sterilization on any moving joints.

### Sterilization and Packaging


- Prior to sterilization, apply a surgical grade instrument milk/anticorrosive lubricant on any moving joints. Do not use without applying anticorrosive lubricant onto moving parts as galling could occur and may cause damage to the instrument.
- Use a validated, properly maintained steam sterilizer. Always follow instructions of the machine manufacturer.
- Do not exceed temperatures over 140°C (284°F) during the sterilization cycle and remove the cleaning wire before sterilizing the devices.
- Instruments and case/tray configurations should be double wrapped or packaged following hospital procedures and as defined by ANSI/AAMI ST 46-1993. The packaging for terminally sterilized reusable instruments should be suitable for steam sterilization and the appropriate grade for the weight of the instruments. Additionally, the blue wrap should be compliant to the following requirements: AAMI ST79, ISO 11607 and FDA 510(k) clearance.
- The process parameters are validated and recommended for sterilization. Steam autoclave (moist heat) sterilization using a pre-vacuum cycle is recommended, as noted above. Autoclaves should comply with the requirements of, and be validated and maintained in accordance with, EN285, EN13060, EN ISO17655 and ANSI/AAMI ST79. And after sterilization, reusable instruments should be stored in the sterilization wrap in a dry and dust-free location.
- The shelf life is dependent on the sterile barrier that is used, storage method, environmental conditions, and handling. A maximum shelf life for sterilized reusable instruments should be defined by each healthcare facility based on the recommendations of the wrap manufacturers.
- Devices shall be stored in a dry, clean environment at an ambient room temperature.
- Effective sterilization can be achieved following the steam cycle:

Sterilization Methods		
Pre-vacuum	Minimum Temperature	Minimum Temperature
	132.2°C/270°F	135°C/275°F
Full Cycle Time	4 Minutes	3 Minutes
Drying Time	30 Minutes	50 Minutes

### Storage and Product Life

- Service of life of this device: 5 years
- Always inspect the devices prior to and after use to confirm proper functioning.
- Do not store this device in high temperature areas or in areas with high humidity where the temperature has dramatic variations.

**Note:** If a used device is returned, it must be decontaminated. Contact Mizuho America via phone (800) 699-CLIP (2547) or email [customerservice@mizuho.com](mailto:customerservice@mizuho.com) to obtain an RGA for return of any damaged instruments.

	Manufacturer		Fragile		Read the Instructions
	Attention		Non Sterile		

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