


## Instruction for Use

Trade Name: Raised armboard

### Warning

1. When using operating table accessories, always pay close attention to the patient. Patients left in the same posture for an extended period of time are at risk for developing bedsores.
2. When using this product for a patient in the lateral recumbent position, secure the arm on the lower side onto the Armboard.  
Failure to secure the arm to the Armboard may lead to a patient injury.




### Contraindication/Prohibition

1. This product is only validated for use with Mizuho operating tables.
2. Ensure clamping handle is tight.  
Failure to do so may result in a patient injury.

### Shape/ Structure

#### 1. Shape/Structure



Code No.	Product Description	A (mm)	Weights (kg)
08-066-06	Raised armboard with clamp (730 mm)	730	3.3
08-066-10	Raised armboard with clamp	500	3.0
* 08-066-10-UE	Raised armboard (UE)	450	2.0
* 08-110-01	Clamp for round bar	—	0.9
* 08-117-09	Clamp for round bar (R1)	—	0.9
* 08-117-11	Clamp for round bar (R2)	—	0.9
08-110-20	Universal clamp for round bar	—	1.0

### 2. Material

- (1) Main Body  
 Armboard : Steel plate, Urethane foam  
 Support pole : Stainless steel rod  
 Arm strap : Nylon resin
- (2) Clamp: Stainless steel


### Intended purpose

This is an accessory to an operating table.  
 This is a tool for supporting the arm on the upper side of a patient lying in the lateral recumbent position for a thoracotomy or similar procedure.

### Specifications

Item	Specification
Raised Armboard	<ul style="list-style-type: none"> <li>· Mounts to the clamp</li> <li>· Is secured using the clamp handle</li> <li>· Is adjusted to desired angle by loosening the clamp handle</li> </ul>
Clamp	<ul style="list-style-type: none"> <li>· Mounts to side rail of an operating table.</li> </ul>

### Instructions for use



- ① Clamp for round bar      ④ Armboard
- ② Clamp handle            ⑤ Armboard ball joint handle
- ③ Support pole            ⑥ Arm strap

1. Attach the clamp for round bar① to the side rail of an operating table.
2. Insert the support pole③ into the clamp for round bar① and tighten the clamp handle② to temporarily secure the support pole.
3. Move the support pole③ to a position appropriate for the patient's posture.
4. After the adjustment, rotate the clamp handle② to secure the clamp for round bar① and the support pole③.
5. Loosen the Armboard ball joint handle⑤ and adjust the angle of the Armboard④.
6. After the adjustment, tighten the Armboard clamp handle⑥ to secure the Armboard④.
7. Secure the patient's arm with the arm strap⑥.

### Caution

1. Important caution  
 When working with an operating table, take care not to allow this product to make contact or interfere with the table top or with other tools and appliances used in combination with this product.

**Storage/Life**

1. Do not store the device in high temperatures or in areas with high humidity where the temperature has drastic variations.
2. Service life of this product: 7 years  
(subject to the specified maintenance and inspection and is stored properly)

**Maintenance/Inspection**

1. Check that this product is not damaged or broken before and after use.
2. Cleaning and disinfection  
Wipe off blood, chemicals, contaminants, and other stains with water, and clean the device with gauze or other material moistened with alcohol to disinfect.
3. In case of a malfunction  
If you need to have this product repaired, clearly label the device as "Broken", "Need repair" etc. and contact your local dealer or Mizuho.

<By suppliers>

Maintenance and inspection can only be carried out certified suppliers of Mizuho.

**Packing**

1 piece per pack

**Matters related to warranty period**

MIZUHO Corporation will repair defective parts of this product without charge for one year from the date of delivery/installment except for cases of damage caused by a third party's repair, act of nature, improper use or damage on purpose. All other warranty terms and conditions are subject to regulations of MIZUHO Corporation

**\*\* Name and address of manufacturer**

MIZUHO Corporation  
3-30-13 Hongo, Bunkyo-ku, Tokyo 113-0033, Japan  
<http://www.mizuho.co.jp>

**\*\* Authorized Representative Europe**

Emergo Europe B.V.  
Prinsessegracht 20, 2514 AP, The Hague  
The Netherlands