MIZUHO

CE
0123

SUGITA

AVM Microclips

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INSTRUCTIONS FOR THE USE OF SUGITA AVM MICROCLIPS

Index

- | | |
|---------------------------------------|---------------------------------------|
| 1. Specifications | 10. Warning |
| 2. Types | 11. Caution |
| 3. Spring Type | 12. Defects and Adverse Events |
| 4. Action Mechanism | 13. Storage - Care and Handling |
| 5. Surface Finish | 14. Sterilization |
| 6. Blade Cross-sectional Shape | - Wrapping Prior to Sterilization - |
| 7. Indication for Use | - Validated Sterilization Procedure - |
| 8. Directions for Use | - Storage - |
| 9. Contraindications and Prohibitions | |

1. Specifications

Sugita AVM Microclips are made of cobalt-chromium-nickel-molybdenum-iron alloy (ISO 5832-7/ASTM F1058 G2). “Sugita AVM Microclip” herein referred to as “clip”.


ASTM: American Society for Testing and Materials

ISO: International Organization for Standardization

This clip designed for **single-use only** is exclusively indicated only for use with the Sugita AVM Microclip Applying Forceps (hereinafter this is called “clip applier”), and requires sterilization prior to use. Reuse may cause the breakage of clips and lead to infection.

The closing force of each clip is measured during the final quality assurance stage with a computerized electronic scale at a measuring point located 1.0 mm from the distal end of the blade length (L) of the clip [**Figure 1**] while its blades are opened at a 1.0 mm width. The result of this measurement is shown in newtons and grams on the package label.

The product information including name, size are as indicated on the package label.

[Figure 1]

2. Types

Prod. Code	Length of blade
07-932-02	2 mm
07-932-03	3 mm
07-932-04	4 mm
07-932-05	5 mm

For specific dimensions and options, please see this device's catalog.

3. Spring Type : Integral Single Bend

4. Action Mechanism : Alligator Action

5. Surface Finish

- Clip Body: Satin
- Material Treatment: Electrolytic Polish

6. Blade Cross-sectional Shape : Rolled

7. Indication for Use

This device is indicated for temporary occlusion of cerebral blood vessels.

For the US market

This device is indicated for temporary occlusion of blood vessels.

8. Directions for Use


This device is packaged non-sterile, and it must be properly sterilized before use.

Each clip must be held with the specified clip applicator.

When holding the clip, it must be verified that the clip is correctly placed in the holding jaws of the clip applicator as shown in **[Figure 2]**.

Improper holding of the spring of clip may cause damage or ejection

from the clip applier. As this may present a risk during surgery, correct holding of the clip must always be kept in mind.

[Figure 2]

- Verify correct placement of the clip in the recesses of the holding jaws of clip applier.
- Do not hold the clip in angled positions.
- If the spring section of the clip is not properly held with the holding jaws of clip applier, the clip may slip or spring out from the jaws. This may cause damage to the closing force of the clip, which will prevent closure of the targeted cerebral blood vessels.

It is recommended hospital personnel records the type, lot and serial numbers of the clip used for each patient.

To support traceability, Patient/Trace labels are included in the clip package.

9. Contraindications and Prohibitions

1) Prohibition of reuse

This device is for single-time use only. Once the blades of clip have been opened or the clip has come into contact with blood or tissue of a patient, it must be discarded. Reuse may cause the breakage of clips and lead to infection. [Cause of breakage/risk of infection]

2) Prohibition of use for unintended purposes

This device should not be used for purposes other than indication for use.

The clips should not be used for permanent occlusion of, or permanent placement in blood vessels.

3) Use only the exclusive clip applier when handling this clip to avoid changes in the clips closing forces or mechanical features of the clip.

Do not use other manufacturers' clip appliers with this clip as there is risk of harm to the patient, breakage or damage to the clip.

4) Do not open the blades of clip directly with fingers. This may cause corrosion or alter the proper holding force. The blades of clip should be opened only when spring of the clip is placed in the proper position inside the clip applier.

5) Modifying or adulteration of this product is not allowed.

6) In order to preclude any possible Galvanic reaction, do not use this device with implants of other manufacturers or with implants made of different materials. These Galvanic reactions may cause corrosion or alter the proper holding force.

-
- 7) This device must not be used in patients with serious allergies to metals or exogenous materials.
 - 8) This device must not be used in patients with infectious diseases or sepsis as a complication.
 - 9) This device must not be sterilized using low-temperature hydrogen peroxide gas plasma.
This product is not compatible with sterilization using low-temperature hydrogen peroxide gas plasma. [Sterilization using such method may cause surface color deterioration and affect its characteristics.]

10. Warning

- 1) It is recommended the patient is informed about the procedure and characteristics of this clip.
- 2) Handling of this product:
Rough handling of this product may cause deformation or damage to the clip and may harm the patient.

11. Caution

Use of this product is restricted to surgeons with proper training in the use of this product.

Impact should not be applied to its package or aluminum container in order to avoid the dropout of clips.

-
- The aluminum container is sealed to maintain the product integrity. Clips removed from their original packaging or sterilized by user cannot be returned.
 - This device can be re-sterilized only if the sterilization bag or wrap is not opened in the operation room.
 - If the sterilization bag or wrap is opened in the operating room, the product must be discarded.
 - Provide a sterile barrier (pouch or wrap) to protect the clip and its container from contamination. Once the sterile barrier has been removed and the clip and its container has been subjected to possible contamination, the clip and its container must be disposed.

MR Conditional

Non-clinical testing demonstrated that this device is MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5 tesla or 3 tesla, only
- Maximum spatial gradient magnetic field of 1,000 gauss/cm (10-T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per

pulse sequence) in the Normal Operating Mode
Under the scan conditions defined, this device is expected to produce a maximum temperature rise of 2.4°C after 15-minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by this device extends approximately 25-mm from this device when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

NOTE: Do not take the clip applicator into the MR environment.

This device contains material which has been determined to be MR Unsafe.

12. Defects and Adverse Events

1) Serious adverse events

- Displacement of the clip (slippage)
- Ejection of the clip from clip applicator
- Breakage of the clip
- Cerebrovascular injury
- Hemorrhage
- Death

2) Other adverse events

- Infection of surgical wound
- General surgical complications
- Delayed clip displacement due to the application of magnetic field
- Occlusion of adjacent blood vessel

-
- Foreign body reaction

13. Storage – Care and Handling

This device is prepared in a non-sterile package and should be stored in its package in a cool, dry, and secure location.

Do not use any clip applicator other than the exclusive clip applicator.

Handling of this device should be minimized to reduce contamination or damage.

14. Sterilization

This device is supplied non-sterile. Therefore, this device must be sterilized by the users (healthcare professionals) prior to surgery using validated sterilization processes in accordance with guidelines stipulated by the government.

Proper sterilization and handling of this clip is important for assuring that the clips meet sterile specifications prior to surgery as an aid in minimizing the occurrence of infection in the surgical patient.

Other procedures are also important in minimizing contamination in an operating room including the use of a surgical unit with laminar flow equipment at the surgeon's discretion.

Sterilization equipment should be validated to assure that it sterilizes specific load configurations properly. The end user is responsible to validate the sterilization process in the individual sterilizer and assure the sterility.

Sterilization equipment should be cleaned and maintained in accordance with standardized operating procedures specified by its manufacturer. Also, sterilize the clip applicator.

- Wrapping Prior to Sterilization -

Do not wash the clip. This clip is not designed to be washed or cleaned by the user.

Remove the sealed aluminum container from the outer box/original packaging and wrap the container with sterilization packaging material (nonwoven fabric or sterilization bag/pouch).

Place in an appropriate sterilization bag/pouch or wrap and sterilize.

Do not remove the clip out of the aluminum container.

To support sterilization, each individual clip is retained in a heat-resistant autoclavable foamed silicone material within an aluminum container.

Do not use low-temperature hydrogen peroxide gas plasma sterilization.

- Validated Sterilization Procedure -

1. Load into autoclave.
2. Set the autoclave with the following parameters.
 - Temperature: 134°C
 - Cycle Type: Pre-vacuum, Exposure time: 3 minutes

Do not sterilize using low-temperature hydrogen peroxide gas plasma. It may cause discoloration and/or deformation.

Device must be sterilized by users in accordance with validated sterilization procedures, such as an autoclave, that are regulated by medical organizations in each country or region.

- Storage -

After sterilization, the clip should be stored and protected in a clean and dry area.

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