MIZUHO

CE
0123

SUGITA

Booster CLIPS

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INSTRUCTIONS FOR THE USE OF SUGITA BOOSTER CLIPS

Index

- | | |
|---|---|
| 1 . Specifications | 10. Warning |
| 2 . Types | 11. Use of MRI |
| 3 . Spring Type | 12. Caution |
| 4 . Action Mechanism | 13. Defects and Adverse Events |
| 5 . Surface Finish | 14. Storage - Care and Handling |
| 6 . Blade Cross-sectional Shape | 15. Validated Cleaning and
Sterilization Procedure |
| 7 . Indication for Use | |
| 8 . Directions for Use | |
| 9 . Contraindications and
Prohibitions | |

INSTRUCTION MANUAL FOR SUGITA BOOSTER CLIPS


1. Specifications

Sugita Booster Clips are used to reinforce the closing force of standard and large Sugita Permanent Aneurysm Clips as shown in [Table 1].

The Sugita Booster Clip should be used only in tandem with a previously applied Sugita Aneurysm Clip. **Never apply directly to the cerebral aneurysm or to any blood vessel.** The Sugita Booster Clip must be placed as close to the spring mechanism of the Sugita Aneurysm Clip as possible and must securely grip the proximal part of the clip blades. [Figure 1]

The Sugita Booster Clips should not be used with clips other than the standard or large Sugita Aneurysm Clips.

SUGITA BOOSTER CLIP


[Figure 1]

Application type	Prod. Code
Standard Type	07-940-01~09·64~66·68·69 07-943-01~03
Large Type	07-940-10~18·20~25·47~50 60~63·70~73

[Table 1]

These Sugita Booster Clips are made of cobaltchromium-nickel-molybdenum-iron alloy (ISO 5832-7/ASTM F1058G2).

Each Sugita Booster Clip is designed for **single-**




time use only with the use of Sugita Booster Clip Applying Forceps, which are exclusive to this particular clip. Reuse may cause the breakage of clips and lead to infection.

ASTM: American Society for Testing and Materials; ISO: International Organization for Standardization.

A lot number is shown on each Sugita Booster Clip in order to ensure Clip traceability.

The closing force of each clip is measured during the final quality assurance stage with a computerized electronic scale at a specified measuring point as shown in **[Figure 2]**. The result of this measurement is shown in newtons and grams on the package label.

[Figure 2]

2. Types

○ Permanent Type

Type	Color
Booster	Blade tip-yellow

3. Spring Type: Integral Coil

4. Action Mechanism: Alligator Action

5. Surface Finish

- Clip Body: Satin
- Material Treatment: Electrolytic Polish

6. Blade Cross-sectional Shape: Rolled

7. Indication for Use

Sugita Booster Clips are used to reinforce the closing force of standard and large Sugita Aneurysm Clips for permanent occlusion (permanent clips)

8. Directions for Use


All Sugita Booster Clips are packaged non-sterile and they must be properly sterilized before use.

Each Sugita Booster Clip must be held with specified Sugita Booster Clip Applying Forceps.

When holding the clip, it is important to verify that the clip is correctly placed in the holding jaws of the


applying forceps as shown in **[Figure 3]**.

Improper holding of the clip may cause damage or ejection from the applying forceps. As this may present a risk during surgery, correct holding of the clip must always be kept in mind.

[Figure 3]

- Verify correct placement of the clip in the recesses of the holding jaws of forceps.
- Do not hold the clip in angled positions.
- If the spring of the clip is not properly held with the holding jaws of forceps, the clip may slip or flick off. This may cause damage to the spring section of

the clip with possible fracture of the spring.

It is very important to verify correct placement of Sugita Aneurysm Clips to the target site. It is recommended that the hospital record the type, lot and serial numbers of the clip used for each patient. Patient/Trace labels are included in the clip package for this purpose.


9. Contraindications and Prohibitions

1) Prohibition of reuse

Sugita Booster Clips are for single-time use only.



Once the blades of clip have been opened or the clip has come into contact with blood or tissue of a patient, it must be discarded. Reuse may cause the breakage of clips and lead to infection. [Cause of breakage/risk of infection]

2) Prohibition of use for unintended purposes


The Sugita Booster Clip should not be used for purposes other than reinforcing the closing force of Sugita Permanent Aneurysm Clips (standard and large types).

Do not use the Sugita Booster Clip with Fenestrated, Long, or Mini Sugita clips.

- 
- 3) Use only the exclusive Sugita Booster Clip Applying Forceps when handling the Sugita Booster Clips to avoid the changes in the closing forces or mechanical features of the clip. Do not use Sugita Booster Clips with clip forceps of other manufacturers.
- 4) Do not open the blades of clip directly with fingers. This may cause corrosion or alter the proper holding force. The blades of clip should be opened only when the spring of clip is placed in a proper position of the forceps.
- 5) Sugita Booster Clips should not be bent or
- 
-

remodeled.

- 6) Do not use Sugita Booster Clips with implants of other manufacturers or with implants made of different materials in order to prevent a possible Galvanic reaction. This may cause corrosion or alter the proper holding force.
 - 7) Sugita Booster Clips must not be used in patients with serious allergies to metals or exogenous materials.
 - 8) Sugita Booster Clips must not be used for patients with infectious diseases or sepsis as a complication.
 - 9) Sugita Booster Clips must not be sterilized using low-temperature hydrogen peroxide gas plasma. This product is not compatible with sterilization using low-temperature hydrogen peroxide gas plasma.
[Sterilization using such method may cause surface color deterioration and affect its characteristics.]
- 

10. Warning

1) Full explanation should be given to each patient about the use of Sugita Booster Clips as well as the benefits/risks of surgery involving its use.

2) Closing status of cerebral aneurysm clips should be verified immediately following surgery.

[The most serious adverse events expected to be reported in conjunction with this booster clip may include displacement of the clip, slip-out of the clip, ejection of the clip from clip forceps, clip fracture, laceration of the aneurysm, cerebrovascular injury, and hemorrhage.]

3) Prohibition of secondary processing on this product
Secondary processing on this product in the form of, for example, incision by applying impact or vibration on it should not be done. [It may cause breakage.]

4) Prohibition of rough handling

Rough handling of this product such as handling that may cause deformation or damage should be avoided. [It will markedly shorten the life of the implant.]


11. Use of MRI


MR Conditional

Non-clinical testing demonstrated that the Sugita Booster Clip is MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5 tesla or 3 tesla, only
- Maximum spatial gradient magnetic field of 1,000 gauss/cm (10-T/m)

•Maximum MR system reported, whole body averaged
specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode

Under the scan conditions defined, the Sugita Booster Clip is expected to produce a maximum temperature rise of 2.4°C after 15-minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the Sugita Booster Clip extends approximately 25-mm from this device when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

NOTE: Do not take Sugita Booster Clip Applying Forceps into the MR environment. This device contains material which has been determined to be MR Unsafe.

12. Caution

Use of the Sugita Booster Clips is restricted to surgeons with proper training and experience in the use of intracranial aneurysm clips.

13. Defects and Adverse Events

1) Serious adverse events

- Displacement of the clip (slip-out)
- Ejection of the clip from clip forceps
- Breakage of the clip
- Rupture of a cerebral aneurysm
- Cerebrovascular injury
- Hemorrhage
- Death

2) Other adverse events

- Infection of surgical wound
 - General surgical complications
 - Delayed clip displacement due to the
-

-
- application of magnetic field
 - Occlusion of adjacent blood vessels
 - Foreign body reaction


14. Storage - Care and Handling

Each Sugita Booster Clip is prepared in a nonsterile package and should be stored in its package in a cool, dry, and secure location.

Do not use any clip applying forceps other than the exclusive Sugita Booster Clip Applying Forceps.


15. Validated Cleaning and Sterilization Procedure

Proper sterilization and handling of the Sugita Booster Clips is important for assuring that the clips meet sterile specifications prior to implantation as an aid in minimizing the occurrence of infection in the surgical patient.

Other procedures are also important in minimizing contamination in an operating room including the use of a surgical unit with laminar flow equipment at the surgeon's discretion.

Sugita Booster Clips are supplied non-sterile.

Therefore, Sugita Booster Clips must be sterilized by the users (healthcare professionals) prior to surgery using validated sterilization processes in accordance with guidelines stipulated by the government.

Sugita Booster Clips and instruments must be cleaned to reduce bioburden before first sterilization and resterilization.

Do not use chlorine-based solution when cleaning this product. This could cause metallic corrosion.

1) Cleaning

Put the clip in the designated Sugita Clip Case and close the lid.

< In case of using Ultrasonic cleaner>

- ① Fully immerse the case in the water (tap water) of the washing tank.
 - ② Ultrasonic cleaning.
Cleaning time: 15 minutes or more
Detergent: Neutral enzyme type detergent (Diluted concentration 1.0%)
 - ③ Rinsing
 - Make sump water with running water (tap water or deionized water), immerse the case completely in water, shake the case 5 times up and down.
 - Renew the water and shake the case 5 times up and down in the water again.
 - Perform a series of rinsing operations three or more times in total.
 - ④ Dry thoroughly.
Recommended drying condition: 110 °C 30 minutes or more
-


<In case of using washer disinfector (automatic washing)>

Stage	Procedure	Temperature [°C]	Time [min]	Quality of water	Remarks
I	Preclean	<25	3	Tap water	
II	Washing	50	10	Tap water	Alkaline detergent (such as cleaning agent for aluminum)-Diluted concentration 0.5%
III	Rinsing	<25	4	Deionized water	Set rinsing process program as "1 time (1 minute) × 4 times = 4 minutes"
IV	Hot water disinfection	–	–	Deionized water	Follow the program of the device.
V	Drying	–	–	–	

2) Inspection

Check products visually one by one and do not use if the following abnormalities are found.

- Damage
- Deformation
- Corruption
- Inappropriate blade position (deviation, clearance, tip unevenness, etc.)

•Remaining dirt

3) Packaging

Put the clip in the designated Sugita Clip Case and close the lid.

Then pack the case containing the product using standard medical packaging material (non-woven fabric or sterilized bag) with protective sterile barrier complying with ISO 11607 for steam sterilization

4) Sterilization

Perform sterilization treatment under sterilization condition "134 °C • 3 minutes or more" with a pre-vacuum type autoclave sterilizer which was confirmed to function properly.

Do not use low temperature hydrogen peroxide gas plasma sterilization. The surface fades and affects its properties.

5) Storage

Do not store the device in high temperatures, or in areas with high humidity and/or where the temperature varies.

Please store the device in clean dust-free and non-splashed area.

- General Requirements -

1. Handling of the Sugita Booster Clips should be minimized to reduce contamination.
 2. Sterilization equipment should be validated to assure that it sterilizes specific load configurations properly. The end user is responsible to validate the sterilization process in the individual sterilizer and assure the sterility.
 3. Sterilization equipment should be cleaned and maintained in accordance with standardized operating procedures specified by its manufacturer.
-

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