**MIZUHO**

**CE**  
0123

**SUGITA**

**Titanium ANEURYSM CLIPS**

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**Emergo Europe B.V.**  
Prinsessegracht 20, 2514 AP, The Hague  
The Netherlands

**MIZUHO Corporation**

3-30-13 Hongo, Bunkyo-ku, Tokyo 113-0033, Japan

TEL(03)3815-3096 FAX(03)3813-5068

<http://www.mizuho.co.jp>

CK07-934-01WW

2020-02 Ver. 15

**INDEX**

English	1	Italiano	65
Deutsche	21	Español	87
Français	43	日本語	109

---

# INSTRUCTIONS FOR THE USE OF SUGITA TITANIUM ANEURYSM CLIPS

## Index

- |  |  |
|--|--|
| 1 . Specifications                     | 11. Use of MRI                           |
| 2 . Types                              | 12. Caution                              |
| 3 . Spring Type                        | 13. Defects and Adverse Events           |
| 4 . Action Mechanism                   | 14. Storage - Care and Handling          |
| 5 . Surface Finish                     | 15. Validated Cleaning and Sterilization |
| 6 . Blade Cross-sectional Shape        | Procedure                                |
| 7 . Indication for Use                 | - Validated Sterilization Procedure -    |
| 8 . Directions for Use                 | - Storage -                              |
| 9 . Contraindications and Prohibitions |  |
| 10. Warning                            |  |

---

## 1. Specifications

Sugita Titanium Aneurysm Clips series offer the permanent type and temporary type.

“Sugita Titanium Aneurysm Clips” herein referred to as “clip”.

This device is made of titanium 6-aluminium 4-vanadium alloy (Ti-6Al-4V ELI; as specified by ISO5832-3/ASTM F136) and unalloyed titanium (specified by ISO 5832-2/ASTM F67 G2).

ASTM: American Society for Testing and Materials

ISO: International Organization for Standardization

This clip is designed for **single-use only**, is exclusively indicated only for use with the Sugita Titanium Aneurysm Clips Applying Forceps (hereinafter this is called "clip applicator"), and requires cleaning and sterilization prior to use. Reuse may cause breakage of clips and may cause injury or infection.


In order to ensure the traceability, a lot number is shown on each clip.

---

---

The closing force of each clip is measured during the final quality assurance stage with a computerized electronic scale at a measuring point located  $1/3 L$  from the distal end of the blade length ( $L$ ) of the clip [Figure 1] while its blades are opened at a 1.0 mm width. The result of this measurement is shown in newtons and grams on the package label.

The product information including name and size are as indicated on the device's package label.

[Figure 1]

---

## 2. Types

### Permanent Type

Type	Color
Standard Type	Whole clip-blue
Mini Type	Blade only-blue
Fenestrated Type	Spring-purple
	Remaining whole body-blue

### Temporary Type

Type	Color
Standard Type	Blade tips-yellow
	Remaining whole body-blue
Mini Type	Blade tips-yellow
	Remaining blade-blue

For specific dimensions and shapes, please see this device's catalog.

---

- 
- 3. Spring Type** : Integral Coil  
**4. Action Mechanism** : Alligator Action  
**5. Surface Finish**  
- Clip Body : Satin  
- Material Treatment : Anodizing Treatment  
**6. Blade Cross-sectional Shape** : Rolled

### **7. Indication for Use**

The permanent type is indicated for the permanent occlusion of the neck of cerebral aneurysms.

The temporary type is indicated for the temporary occlusion of the neck of cerebral aneurysms.

For the US market

The permanent type is indicated for the permanent occlusion of the neck of cerebral aneurysms and peripheral blood vessels.

The temporary type is indicated for the temporary occlusion of the neck of cerebral aneurysms and peripheral blood vessels.

---

---


## 8. Directions for Use

**This device is packaged non-sterile, and must be properly sterilized before use.**

Each clip must be held with the specified clip applicator. Each clip applicator is etched to read; e.g. “FOR (TYPE OF CLIP) Only”.

When holding the clip, it must be verified that the clip is correctly placed in the holding jaws of the clip applicator as shown in **[Figure 2]**.

Improper holding of the spring of clip may cause damage or ejection from the clip applicator. As this may present a risk during surgery, ensure correct holding of the clip.

[Figure 2]

- 
- Verify correct placement of the clip in the recesses of the holding jaws of clip applicator.
  - Do not hold the clip in angled positions.
  - If the spring section of the clip is not properly held with the holding jaws of clip applicator, the clip may slip or spring out from the jaws. This may cause damage to the closing force of the clip, which will prevent closure of the targeted cerebral blood vessels.

It is recommended hospital personnel records the type, lot and serial numbers of the clip used for each patient. To support traceability, Patient/Trace labels are included in the clip package.

## **9. Contraindications and Prohibitions**

### 1) Prohibition of reuse

This device is for single-time use only. Once the blades of clip have been opened or the clip has come into contact with blood or tissue of a patient, it must be discarded. Reuse may cause breakage of

- 
- clips and may cause injury or infection.
- 2) Prohibition of use for unintended purposes  
This device should not be used for purposes other than indication for use.
  - 3) Use only the exclusive clip applier when handling this clip to avoid changes in the clips closing forces or mechanical features of the clip. Do not use other manufacturers' clip appliers with this clip as there is risk of harm to the patient, breakage or damage to the clip.
  - 4) Do not open the blades of clip directly with fingers. This may cause corrosion or alter the proper holding force. The blades of clip should be opened only when spring of the clip is placed in the proper position inside the clip applier.
  - 5) Modifying or adulteration of this product is not allowed.
  - 6) In order to preclude any possible Galvanic reaction, do not use this device with implants of other manufacturers or with implants made of different
-

---

materials. These Galvanic reactions may cause corrosion or alter the proper holding force.

- 7) This device must not be used in patients with serious allergies to metals or exogenous materials.
- 8) This device must not be used in patients with infectious diseases or sepsis as a complication.
- 9) This device must not be sterilized using low-temperature hydrogen peroxide gas plasma.

This product is not compatible with sterilization using low-temperature hydrogen peroxide gas plasma. [Sterilization using such method may cause surface color deterioration and affect its characteristics.]

## **10. Warning**

- 1) It is recommended the patient is informed about the procedure and characteristics of this clip.
- 2) Closing status of the blades of clip should be verified immediately following surgery.  
[The most serious adverse events expected to be reported in conjunction with the aneurysm clip may

---

include displacement of the clip, slip-out of the clip, ejection of the clip from clip applier, clip fracture, laceration of the aneurysm, cerebrovascular injury, and hemorrhage. In addition, in consideration of the structure of the clip, risk of failures/adverse events may be increased by applying only the tips of clips, especially the ones with long blades, to the target tissue.]

3) Handling of this product:

Rough handling of this product may cause deformation or damage to the clip and may harm the patient.

---

## 11. Use of MRI

### MR Conditional

Non-clinical testing demonstrated that this device is MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5 tesla or 3 tesla, only
- Maximum spatial gradient magnetic field of 1,000 gauss/cm (10-T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode

Under the scan conditions defined, this device is expected to produce a maximum temperature rise of 2.2°C after 15-minutes of continuous scanning (i.e.,

---

---

per pulse sequence).

In non-clinical testing, the image artifact caused by this device extends approximately 5-mm from this device when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

**NOTE:** Do not take the clip applicator into the MR environment.

The clip applicator contains material which has been determined to be MR Unsafe.

## **12. Caution**

Use of this product is restricted to surgeons with proper training in the use of this product.

## **13. Defects and Adverse Events**

1) Serious adverse events

- Displacement of the clip (slippage)
- Ejection of the clip from clip applicator
- Breakage of the clip
- Rupture of a cerebral aneurysm

- 
- Cerebrovascular injury
  - Hemorrhage
  - Death
- 2) Other adverse events
- Infection of surgical wound
  - General surgical complications
  - Delayed clip displacement due to the application of magnetic field
  - Occlusion of adjacent blood vessel
  - Foreign body reaction

## **14. Storage – Care and Handling**

This device is prepared in a non-sterile package and should be stored in its package in a cool, dry, and secure location.

Do not use any clip applicator other than the exclusive clip applicator.

Handling of this device should be minimized to reduce contamination or damage.

---

## **15. Validated Cleaning and Sterilization Procedure**

This device is supplied non-sterile. Therefore, this device must be sterilized by the users (healthcare professionals) prior to surgery using validated sterilization processes in accordance with guidelines stipulated by the government.

Proper sterilization and handling of this clip is important for assuring that the clips meet sterile specifications prior to surgery as an aid in minimizing the occurrence of infection in the surgical patient.

Other procedures are also important in minimizing contamination in an operating room including the use of a surgical unit with laminar flow equipment at the surgeon's discretion.

Sterilization equipment should be validated to assure that it sterilizes specific load configurations properly. The end user is responsible to validate the sterilization process in the individual sterilizer and assure the sterility.

---

---

Sterilization equipment should be cleaned and maintained in accordance with standardized operating procedures specified by its manufacturer.

Clean and sterilize the clip applicator according to the clip applicator instruction manual.

This clip and instruments must be cleaned to reduce bioburden before first sterilization and re-sterilization.

Do not use chlorine-based solutions when cleaning this product. This could cause metallic corrosion.

#### 1) Cleaning

Put the clip in the designated Sugita Clip Case or Sugita Clip Tray (hereinafter this is called "clip case or tray" ) and close the lid.

<In case of using Ultrasonic cleaner>

(1) Fully immerse the case (tray) in the water (tap water) of the washing tank.

(2) Ultrasonic cleaning

Cleaning time: 15 minutes or more

Detergent: Neutral enzyme type detergent (Diluted)

---

concentration 1.0%)

(3) Rinsing

- Make sump water with running water, immerse the case (tray) completely in water (tap water or deionized water), shake the case (tray) 5 times up and down.
- Renew the water and shake the case (tray) 5 times up and down in the water again.
- Perform a series of rinsing operations three or more times in total.

(4) Dry thoroughly.

Recommended drying condition: 110 °C, 30 minutes or more

---

<In case of using washer disinfectant (automatic washing)>

Stage	Procedure	Temperature [°C]	Time [min]	Quality of water	Remarks
I	Preclean	<25	3	Tap water	
II	Washing	55	10	Tap water	Alkaline detergent (such as cleaning agent for aluminum) that does not affect anodic oxidation coatings Diluted concentration 0.5%
III	Rinsing	<25	4	Deionized water	Set rinsing process program as "1 time (1 minute) × 4 times = 4 minutes"
IV	Hot water disinfection	—	—	Deionized water	Follow the program of the device.
V	Drying	—	—	—	

## 2) Inspection

Check products visually one by one and do not use if the following abnormalities are found.

- Damage
- Deformation
- Corruption
- Inappropriate blade position (deviation, clearance, tip unevenness, etc.)
- Remaining dirt

---

Since the color of the surface may be faded in repetitive processing, do not use products that cannot be clearly distinguished with color coding.

### 3) Packaging

Place clip in the designated clip case or tray and close the lid.

Then pack the case (tray) containing the product using standard medical packaging material (non-woven fabric or sterilized bag) with protective sterile barrier complying with ISO 11607 for steam sterilization.

### **- Validated Sterilization Procedure -**

With the pre-vacuum type high-pressure steam sterilizer which functions properly, this product shall be sterilized so that sterilization conditions [134°C for 3 minutes or more] are maintained.

Do not sterilize using low-temperature hydrogen peroxide gas plasma. It may cause discoloration and/or deformation.

---

With a high-pressure steam sterilizer that has been confirmed to function properly, the validation of sterilization is verified by a medical institution, sterilization conditions with a proven effectiveness shall be applied.

**- Storage -**

Do not store the device in high temperatures, or in areas with high humidity and/or where the temperature varies. Please store the device in clean dust-free and dry area.

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