MIZUHO

SUGITA

Titanium ANEURYSM CLIPS II BOOSTER

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2014.1(2)

INDEX

English	1	Italian	63
German	21	Spanish	83
French	43	Japanese	105

DIRECTION FOR USE OF THE SUGITA TITANIUM CLIP II BOOSTER

Table of Contents:

- | | |
|---------------------------------------|-------------------------------------|
| 1. Specifications | 10. Warnings |
| 2. Types | 11. Caution |
| 3. Spring Type | 12. Failures and Adverse Events |
| 4. Action Mechanism | 13. Storage Care and Handling |
| 5. Surface Finish | 14. Sterilization Procedure |
| 6. Cross-Sectional Blade Shape | - Cleaning Prior to Sterilization - |
| 7. Indications for Use | - Wrapping Prior to Sterilization - |
| 8. Directions for Use | - Steps - |
| 9. Contraindications and Prohibitions | - General Requirements - |


Instruction Manual for the Sugita Titanium Clip II Booster

1. Specifications:

The Sugita Titanium Clip II Booster increases closing force when bound together with a permanent-type, Standard Sugita Titanium Clip II, as listed in [Table 1].

The Sugita Titanium Clip II Booster has a channel engraved into the tip of the inner blade. Securely attach this channel onto the base of the blade of the Standard Sugita Titanium Clip II whose closing force you want to increase. The two joined clips must be positioned with parallel springs [Fig. 1 (A)]. Do not use if the springs are not in a parallel position [Fig. 1 (B)].

Do not use the Sugita Titanium Clip II Booster with any clip other than the Standard Sugita Titanium Clip II.

[Fig. 1]

[Table 1]

	Product Code
Applicable clips	17-001-** 01, 02, 05, 06, 09 to 11, 15, 17, 18, 20 to 25, 47 to 50, 60 to 63, 65, 66, 69, 70, 73 17-002-01 to 03

The Sugita Titanium Clip II Booster is made of Titanium alloy material, Ti-6Al-4V ELI/ASTM F136 (ISO5832-3), which features metallurgical and biocompatibility properties.


The Sugita Titanium Clip II Booster can be used with the applier for the Standard Sugita Titanium Clip II, and are for **single use only**.

Sugita Titanium Clip II Booster is manufactured

using non-ferromagnetic materials (ASTM F136/ISO 5832-3). Having been verified through tests defined by the ASTM (American Society for Testing and Materials) standards, the clip can be used clinically with 3T (Tesla) or less MRI devices without any adverse effects.

Each Sugita Titanium Clip II Booster is labeled with a manufacturing serial number to ensure traceability.

In the final quality assurance stage, the closing force of clip is measured with a computerized electronic scale at the specified location shown in **[Fig. 2]**. The results of this measurement are shown on the package label in newtons and grams.

[Fig. 2]

2. Types:

○ Permanent-Type

Type	Color
Booster	Spring - Green
	Blade - Green

3. Spring Type: Integral Coil

4. Action Mechanism: Alligator Action

5. Surface Finish:

- Clip Body: Satin

- Material Treatment: Anodizing Treatment

6. Cross-Sectional Blade Shape: Rolled

7. Indications for Use:


The Sugita Titanium Clip II Booster is used to increase closing force when combined with a permanent-type, Standard Sugita Titanium Clip II.

8. Directions for Use:

The Sugita Titanium Clip II Booster is packaged non-sterile; therefore, it must be properly sterilized before use.

The Sugita Titanium Clip II Booster must be applied with the applicator for Standard Sugita Titanium Clip II.

When holding the clip with the applying forceps, verify that the clip is correctly placed in the holding jaws, as shown in [Fig. 3]. Improper hold of the clip may cause damage, slip off or release from the applying forceps. As this may present a risk in surgery, always hold the clip correctly.

[Fig. 3]

-
- Verify the correct placement of the clip in the forceps' holding jaws.
 - Do not hold the clip in any angled position.
 - If the spring of the clip is not properly held within the holding jaws, the clip may slip or flick off. This may cause damage to the spring of the clip and possibly fracture the spring.

It is extremely important to verify that the clip is correctly attached to the cervical segment of the cerebral aneurysm and peripheral blood vessel. It is recommended that surgeons record the clip type used and the manufacturing serial number in the operative notes of the patient. The clip case contains labels for recording purposes.

9. Contraindications and Prohibitions:

1) Do Not Re-Use.

The Sugita Titanium Clip II Booster should be used a single time only.

Once grasped with the clip applying forceps or having come in contact with the patient's blood or tissue, the clip must be discarded, due to possible breakage and a risk of infection.

2) Prohibitions Other than Indications for Use

The Sugita Titanium Clip II Boosters should not be used for any purpose other than to increase closing force in combination with the permanent-type, Standard Sugita Titanium Clip II.

Nor should it be used with clips other than the Standard Sugita Titanium Clip II that are listed in **[Table 1]**.

-
- 3) Exclusive clip applying forceps should be used to ensure that the closing force and mechanical properties of the clip are not altered. Do not use the clip applying forceps of other manufacturers.
 - 4) The clip should not be opened by fingers directly. The clip should only be opened when placed properly within an appropriate clip applying forceps.
 - 5) The clip should not be bent or altered.
 - 6) The clip should not be used with other manufacturers' implants or with an implant made of a different material.
 - 7) Do not use in patients that have severe allergies to metals or allergies to foreign bodies.
 - 8) Do not use in patients that have infections or have a risk of infection.
-

-
- 9) Do not use low temperature hydrogen peroxide gas plasma sterilization.

This product is not suitable for low temperature hydrogen peroxide gas plasma sterilization. (It may cause the surface to fade or may effect its characteristics.)

10. Warnings:

- 1) The use of the Sugita Titanium Clip II Booster, and the benefits and risks arising from its use in surgery should be extensively explained to each patient prior to surgery.
- 2) Immediately after surgery, verify that the cerebral aneurysm clip is properly closed.

Major adverse events that can be expected with the use of this product include: clip migration or slip-

out; dropped clips from clip applying forceps; clip damage; rupture of a cerebral aneurysm; injury to cerebral blood vessels; and bleeding. Additionally, the likelihood of the onset of failures and adverse events increases due to the clip structure, especially when clipping with only the tip of a long-type blade.

3) Secondary processing of the clip is prohibited.

No secondary processing to the surface of the clip, such as etching a mark by striking or using oscillation should be performed. (This may cause breakage.)

4) Rough handling is prohibited.

Do not handle the clip roughly; including, but not limited to altering its shape or scratching it. (This will significantly decrease the life of the implant.)

11. Caution:

Only surgeons with proper training and experience in the use of aneurysm clips should use the Sugita Titanium Clip II Booster.

12. Failures and Adverse Events:

1) Major Adverse Events

- Clip Detachment (Slip-out)
- Clip dropout from applying forceps
- Clip damage
- Rupture of cerebral aneurysm
- Injury to cerebral blood vessels
- Bleeding
- Sudden death

2) Other Adverse Events

- Infection of the surgical wound

-
- General surgical complications
 - Delayed clip detachment due to the application of a magnetic field
 - Occlusion of adjacent blood vessels
 - Reaction to a foreign body

13. Storage Care and Handling:

Each Sugita Titanium Clip II Booster is supplied in a non-sterile package and should be stored in a cool, dry, and secure location.

When handling the clip, use only the applicator for Standard Sugita Titanium Clip II.

Authorized EC-representative:

Emergo Europe

Molenstraat 15 2513 BH The Hague

The Netherlands

Tel : (31)(0)70 345-8570

Fax : (31)(0)70 346-7299

Areas other than EC:

Please contact your local distributors.

14. Sterilization Procedure:

Proper sterilization and handling of the Sugita Titanium Clip II Booster is important for ensuring that the clip meets sterile specifications prior to implantation as an aid in minimizing the occurrence of infection in surgical patients.

Other procedures also important to minimize contamination include the use of laminar flow operating rooms, at the surgeon's discretion.

Sugita Titanium Clip II Booster is supplied non-

sterile. Therefore, they must be sterilized by medical personnel using validated sterilization processes in accordance with guidelines determined by each individual country.

- Cleaning Prior to Sterilization -

Clean the Sugita Titanium Clip II Booster and instruments to reduce the microbial load. Do not use any chlorine-containing solutions as these may cause corrosion.

- Wrapping Prior to Sterilization -

Prior to sterilization, take out the clip from the polyurethane cushion in the clip case, and sterilize the clip only.

An exclusive sterilizing case is available for use.

For high-pressure steam sterilization, the clips should be double wrapped in material for sterile medical packaging.

Also sterilize the applicator for Standard Sugita Titanium Clip II that will be used.

- Steps -

1. Load the double-wrapped sterilization case that contains the Sugita Titanium Clip II Booster into the high-pressure steam sterilizer.
2. Set the temperature to 134°C (sterilization time is 3 minutes or longer) and activate the device.

- General Requirements -

1. After sterilization, the Sugita Titanium Clip II Booster should be stored in a clean place until

-
- used.
2. When wrapping is used, the outer wrapping of the sterilized clip case must be removed after it is placed in the operating room. The inner wrap should remain sealed until just prior to the time of use.
 3. Handling of the Sugita Titanium Clip II Booster should be minimal to reduce the microbial load to the greatest possible extent.
 4. Sterilization equipment should be validated to ensure that it sterilizes specific load configuration properly.
 5. Sterilization equipment should be cleaned and maintained in accordance with standardized operating procedures specified by the sterilizer manufacturer.
-

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