Intended Use

This device is intended for use by a surgeon in a clinical environment for the surgical purpose of grasping and cutting.

CAUTION

Federal law restricts this device to sale by or on order of a physician.

This IFU applies to multiple types of Single Shaft Butterfly Grip Handle Instruments. Prior to use of this product ensure the product is operated and used only by persons with the required training, knowledge and/or experience. No additional training is required as Surgeons have an intimate knowledge of anatomy and physiology, and are trained in the use of similar devices as part of their surgical practice. Read, follow and retain these instructions and use the product only for its intended purpose.

Instruction for Use

Before and after using this device, inspect for function and/or damage, and wash and sterilize in accordance with these instructions. If any issues are found related to the device, do not use the device and contact Mizuho America Representative.

Warnings

- There is a risk of injury and/or failure if the product is not inspected prior to use and after use.
- Risk of injury when applying the product outside your field of view.
- Product should be used according to these instructions and with close attention to patient safety. Failure to follow these instructions and/or adhere to the healthcare facility's requirements may cause harm to the patient and/or healthcare workers.
- Do not reuse the device when it is used on a patient with, or suspected of having, Creutzfeldt-Jakob Disease (CJD) or variant CJD (vCJD).
- The device must be inspected, cleaned, and sterilized before each surgical procedure.

Contraindication/Precautions


- The improper use of this device during the handling, surgical use or reprocessing procedures for which it is indicated, may result in damage or device failure.
- Carefully remove and replace tip protectors as needed.
- Aluminum handles may lose their color over time through normal use and reprocessing.
- To avoid scratching and removal of the ebonized coating, keep ebonized devices separate from other devices.
- Do not use in such a manner as to cause excessive stress or strain at joints; misuse will result in misalignment or cracks at joints or jaws.
- Washing with an improper detergent could result in discoloration or corrosion.
- This device is provided non-sterile.
- Material: Stainless Steel, Aluminum and Titanium.

Maintenance, Inspection and Functional Testing

- Prior to and after each use, and after each cleaning, visually inspect all surfaces for damage, cracks or corrosion and ensure the instrument is functional
- Visually inspect and check all devices for damage and wear; that instrument jaws are free of nicks and damage, and present a continuous edge; that cutting edges are free of nicks and present a continuous edge; that teeth and jaws align correctly; that articulated instruments have a smooth movement without excess play; and that slender devices are not distorted.
- Do not use devices with joints without applying anticorrosive lubricant onto its moving parts as galling could occur and may cause damage to the instrument.
- Do not use if device is defective or damaged.

Decontamination

- Reprocess all devices as soon as they are ready for cleaning; delaying the time for cleaning may damage or stain the device.
- Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion.
- Residues containing chlorine or chlorides (e.g. from surgical residues, drugs, saline solutions, and water for cleaning, disinfection, and sterilization) may cause corrosion and/or surface damage making the product(s) unusable or unrepairable. The product must be removed by rinsing thoroughly with demineralized water and then drying.
- Delicate surgical devices require special handling to prevent damage to the tips. Use caution during cleaning and sterilization. A non-fibrous sponge should be used to wipe off all blood and debris.
- Remove gross contaminants with a steady stream of lukewarm/cool water, not to exceed 35°C (95°F). Rinse each device thoroughly. Do not use saline or chlorinated solutions. Only use distilled or deionized water.
- Do not soak devices in hot water, alcohol, disinfectants or antiseptic solutions. Do not exceed 2 hours soaking in any solution.
- Do not use steel wool, wire brushes or abrasive detergents.
- This product may have a removable cleaning port. The cleaning port is attached to clean the lumen. When cleaning is complete, remove it [Figure 1].
- Devices manufactured from different metals should be processed separately to avoid electrolytic action between the metals resulting in pitting and rusting of devices.
- Open jaws of devices for cleaning. Give special attention to joints, hinges, teeth and serrations. Devices with cleaning ports must be rinsed and cleaned thoroughly.
- The product must be removed by rinsing thoroughly with demineralized water and then drying.
- Follow instructions and warnings as issued by manufacturers of any decontaminants, disinfectants, and cleaning agents used. Wherever possible avoid use of mineral acids or harsh, abrasive agents.
- Use only medical detergents to clean this device.
- Do not use high acid (pH 4.0 or lower) or high alkaline (pH 10 or higher) products for disinfection. Neutral pH detergents 7.0-9.0 are preferred.
- After washing, rinse each device thoroughly. Do not use saline or chlorinated solutions. Only use distilled or deionized water.

Cleaning & Drying: Automated

- Before re-attaching the cleaning port to the product, rinse with tap water and clean with a soft nylon brush until all visible residue is removed from the surface of the instrument. Be especially careful to clean the connecting part of the cleaning port and the product. Upon completion, re-attach the cleaning port to the product.
- Use only validated washer-disinfecter machines and low-foaming, non-ionizing cleaning agents and detergents.
- Follow the manufacturer's instructions for use, which includes warnings, concentrations, and recommended cycles, and load the devices carefully.
- Place heavy instruments on the bottom containers, taking care not to place on delicate devices or overload wash baskets.
- Place instruments with curved surfaces facing down to prevent pooling of water.
- Ensure that soft, high purified water that is controlled for bacterial endotoxins is used in the final rinse stage.
- Devices must be thoroughly dried and all residual moisture must be removed before they are sterilized. Use a soft absorbent towel or cloth to dry external surfaces.
- After washing, apply a water-based, anticorrosive lubricant prior to sterilization on any moving joints.

Sterilization and Packaging

- Prior to sterilization, apply a water-based, anticorrosive lubricant on any moving joints. Do not use without applying anticorrosive lubricant onto moving parts as galling could occur and may cause damage to the instrument.
- Use a validated, properly maintained steam sterilizer. Always follow instructions of the machine manufacturer.
- Do not exceed 140°C (284°F) during sterilization cycle.
- Instruments and case/tray configurations should be double wrapped according to hospital and ANSI/AAMI ST 46-1993 requirements. The packaging for terminally sterilized reusable instruments should be suitable for steam sterilization and the appropriate grade for the weight of the instruments. Additionally, the blue wrap should be compliant to the following requirements: AAMI ST79, ISO 11607 and FDA 510(k) clearance.
- The process parameters are validated and recommended for sterilization. Steam autoclave (moist heat) sterilization using a pre-vacuum cycle is recommended, as noted above. Autoclaves should comply with the requirements of, and be validated and maintained in accordance with, EN285, EN13060, EN ISO17655 and ANSI/AAMI ST79. And after sterilization, reusable instruments should be stored in the sterilization wrap in a dry and dust-free location.
- The shelf life is dependent on the sterile barrier that is used, storage method, environmental conditions, and handling. A maximum shelf life for sterilized reusable instruments should be defined by each healthcare facility based on the recommendations of the wrap manufacturers.
- Devices shall be stored in a dry, clean environment at an ambient room temperature.
- Effective sterilization can be achieved following the steam cycle:

Sterilization Methods		
	Minimum Temperature	Minimum Temperature
Pre-vacuum	132.2°C/270°F	135°C/275°F
Full Cycle Time	4 Minutes	3 Minutes
Drying Time	30 Minutes	50 Minutes

Storage and Product Life

- Service of life of this device: 5 years
- Always inspect the devices between uses to confirm proper functioning.
- Do not store this device in high temperatures or in areas with high humidity where the temperature has dramatic variations.

Note: If a used device is returned, it must be decontaminated. Contact Mizuho America via phone (800) 699-CLIP (2547) or email customerservice@mizuho.com to obtain an RGA for return of any damaged instruments.

Manufacturer
Attention
Read the Instructions
Fragile
Non Sterile

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