



Dr. Youssef's Cranial Nerve Dissection Set

DR. YOUSSEF'S CRANIAL NERVE DISSECTION SET


Comprehensive, high quality, specialized set for nerve dissection

Comprehensive Set for Tumor Treatment



Spatula Dissectors

1.8mm, 2.2mm, 2.4mm, tip sizes



Round Dissectors

2mm size tip with 30°, 45° angles



Putter Dissectors

60°, left and right angles




Recommended intended use for:
acoustic neuroma, meningiomas, other schwannomas.



Teardrop Probes

15°, 45°, 90° angle with ball tip



Needle Dissectors

Straight, 45°, and 90° angles



Cup Curettes

45° & 90° angles, 2x2.5mm & 1.5x2mm cup sizes



Teardrop Back-cutting Knife

90° angle with ball tip



Dr. Youssef's Cranial Nerve Dissection Set

SUPERIOR DESIGN

Comprehensive set for skull base and ENT tumor treatment

Lower profile instruments to accommodate minimally invasive procedures

Extended working length to adapt to a full range of skull base procedures



Uniquely designed tips to mitigate cranial nerve trauma

Vividly colored handles for better instrument identification

Gentle manipulation of tumors

ORDERING INFORMATION

Product #	Description
10000MA	Youssef Cranial Nerve Dissection Set
10001MA	Spatula Dissector, 10cm WL, 2.4mm Tip
10002MA	Spatula Dissector, 10cm WL, 2.2mm Tip
10003MA	Spatula Dissector, 10cm WL, 1.8mm Tip
10004MA	Round Dissector, 10cm WL, 30 Degree Up, 2.0mm Tip
10005MA	Round Dissector, 10cm WL, 45 Degree Up, 2.0mm Tip
10006MA	Putter Dissector, 10cm WL, 60 Degree Up, Right Angle
10007MA	Putter Dissector, 10cm WL, 60 Degree Up, Left Angle
10008MA	Teardrop Probe, 10cm WL, 15 Degree Up, Ball tip
10009MA	Teardrop Probe, 10cm WL, 45 Degree Up, Ball tip
10010MA	Teardrop Probe, 10cm WL, 90 Degree Up, Ball tip
10011MA	Needle Dissector, 10cm WL, Straight
10012MA	Needle Dissector, 10cm WL, 45 Degree Up
10013MA	Needle Dissector, 10cm WL, 90 Degree Up
10014MA	Teardrop Backcutting Knife, 10cm WL, 90 Degree Up, Ball tip
10015MA	Cup Curette, 10cm WL, 45 Degree Up, 2x2.5mm Cup size
10016MA	Cup Curette, 10cm WL, 90 Degree Up, 2x2.5mm Cup size
10017MA	Cup Curette, 10cm WL, 45 Degree Up, 1.5x2mm Cup size
10018MA	Cup Curette, 10cm WL, 90 Degree Up, 1.5x2mm Cup size
10019MA	Cup Curette, 10cm WL, 100 Degree Reverse Angle, 2x2 5mm Cup Size
10020MA	Cup Curette, 10cm WL, 100 Degree Reverse Angle, 1.5x2mm Cup Size



NOTES
