



## Evans 360° Rotatables

Endoscopic and Cranial  
Instrumentation

# Evans 360° Rotatables

Endoscopic and Cranial Instrumentation

## Instrument Set Benefits

- Ideal for endoscopic approaches to pituitary and micro neurosurgical procedures
- 11cm & 15cm lengths available based on surgeon preference
- Newly added detachable flush port facilitates improved cleaning

## Micro Grasping Forceps



*Angled Up*



*Straight*

## Micro Needle Holders



*Straight*



*Curved*



*360° Fully Rotatable Shaft*



## Superior Design

- 360° fully rotatable shafts
- Axial bayonet for improved visualization
- Round handles for easier grip
- Balanced for comfort
- Safely prevents rotation when jaws are closed

## Micro Scissors

---



*Straight*



*Curved*



*125° Back Cutting*



*45° Up*



*Curved, Extended Blades*

## Micro Biopsy Cup Forceps

---



*2x5mm Oval Cup, Straight*



*2x5mm Oval Cup, 45° Up*



*2mm Round Cup, Straight*



*2mm Round Cup, 45° Up*




*4mm Round Cup, Straight*



*4mm Round Cup, Angled Up*

### Ordering Information



#### MICRO NEEDLE HOLDERS

11cm WL	15cm WL	Specifications
07-829-01	07-829-02	Straight
07-829-03	07-829-04	Curved

#### MICRO BIOPSY CUP FORCEPS

11cm WL	15cm WL	Specifications
07-829-11	07-829-12	2mm Round Cup, Straight
07-829-13	07-829-14	2mm Round Cup, 45° Up
07-829-15	07-829-16	2x5mm Oval Cup, 45° Up
07-829-17	07-829-18	2x5mm Oval Cup, Straight
07-829-19	07-829-20	4mm Round Cup, Straight
07-829-21	07-829-22	4mm Round Cup, Ang Up



#### MICRO GRASPING FORCEPS

11cm WL	15cm WL	Specifications
07-829-31	07-829-32	Straight
07-829-33	07-829-34	Angled Up



#### MICRO SCISSORS

11cm WL	15cm WL	Specifications
07-829-41	07-829-42	45° Up
07-829-43	07-829-44	Straight
07-829-45	07-829-46	Curved
07-829-47	07-829-48	125° Back Cutting
NEW 07-829-49MA	07-829-50MA	Curved, Extended Blades



#### TRAY

02-082-00MA	Evans Rotatable Instruments Tray
-------------	----------------------------------

#### FLUSH PORT

P07-829-00	Detachable Flush Port
------------	-----------------------

### Notes

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---