



Kelly Endonasal Instrument Set

All-in-One Endonasal
Transsphenoidal

Designed to support surgical needs for
MICROSCOPIC APPROACHES
WITH ENDOSCOPIC ASSISTANCE

- Low-profile devices that maximizes visibility at the surgical site
- Comprehensive set saves time with less hassle
- Designed in collaboration with Dr. Daniel F. Kelly, M.D.

Product Line-up

- Cottle Nasal Speculum
- Papavero Speculum
- Papavero Crank Handle
- Bayonetted Chisel
- Mallet
- Kerrison Rongeurs
- Mucosal Dissectors
- Suctions
- Decker Cup Forceps
- Yasargil Forceps
- Pistol Grip Micro Scissors
- Ring curettes
- Hemostat
- Weitlaner Retractor
- Mayo scissors
- Metzenbaum scissors
- Spoon Forceps
- Feather Blades and Handle
- Kelly Endonasal Probe, 8Mhz
- Titanium Speculum (optional)




KELLY ENDONASAL INSTRUMENT SET

TRAY 1 - OPENING




Cottle Nasal Speculum

- Short & medium
- 76x7mm & 90x7mm



Papavero Speculum

- Pediatric, Short, Medium & Long
- 65x10mm, 70x13mm, 80x13mm, & 90x13mm




Bayonetted Chisel & Mallet

- 24cm Working Length
- 3mm Jaw
- Mallet with Interchangeable Head



Kerrison Rongeurs

- 180mm Working Length
- 1mm, 2mm & 3mm Up Biting
- 40° Up



Mucosal Dissectors


- Available in Left & Right
- Spear-like Tip



One-Piece Suctions

- 5Fr, 6Fr, 7Fr, 8Fr, 9Fr & 10Fr
- Long (150mm) & Extra Long (170mm) Working Lengths

07-200-01MA
 Kelly Endonasal Set
 Organization Tray #1



TRAY 2 - OPERATING & CLOSING



Spatula Dissectors

- 2.4mm Tip
- Curved Up & Down



Decker Cup Forceps

- Oval & Round Cups
- 2mm, 2x5mm & 2x6mm Cup Sizes
- Straight & 45° Up



Yasargil Forceps

- Bayonet
- 23.5cm Total Length



Adson Tissue Forceps

- 12cm Total Length
- 1x2mm Teeth



Dressing Forceps

- 16cm Total Length



Meningioma Scissors

- 9cm & 11cm Working Length
- Curved & Straight Tips



Pistol Grip Micro Scissors

- 14.5cm Working Length
- Straight & Curved Left & Right



Spoon Forceps

- Cushing, Love-Gruenwald, Spurling Types
- 178mm Working Length
- 2x10mm, 3x10mm & 4x10mm Cup Sizes
- Straight Angle



Hemostat


- 15cm Total Length



Dissectors

Probe teardrop

- 2.5mm Tip
- Straight, 45°, 60° & 90°
- Straight, Up & Down Angles



Ring Curettes

- 1.5mm, 3mm, 5mm & 7mm Cup Sizes
- 45° & 90° Up
- 45° Down



Mezenbaum Scissors

- 18cm Total Length
- Straight & Curved




Mayo Scissors

- 17cm Total Length
- Straight



Weitlaner Retractor

- 13cm Total Length



07-200-02MA
Kelly Endonasal Set
Organization Tray #2

MIZUHO

