



## Kelly Endoscopic Endonasal Instrument Set

# Designed to support surgical needs for ENDOSCOPIC APPROACHES

- Low-profile devices that maximize visibility at the surgical site
- Comprehensive set saves time with less hassle
- Designed in collaboration with Dr. Daniel F. Kelly, M.D.

## Product Line-up

- Cottle Elevator
- Probes
- Nasal Forceps
- Antrum Punch
- Nasal Scissors
- Suctions
- Antrum Cannula
- Weil Blakesley Ronguers
- Kerrison Rongeurs
- Rhoton Dissectors
- Decker Cup Forceps
- Pistol Grip Micro Scissors
- Spoon Forceps
- Ring Curettes
- Hemostat
- Mayo Scissors
- Weitlaner Retractor
- Metzenbaum Scissors
- Spoon Forceps
- Feather Blades and Handle

## Tray 1



### Cottle Elevators

- Double-ended
- 20cm Total Length




### Probes

- Double-ended
- 19cm Total Length



### Rhoton #3 Dissectors

- 12cm Working Length
- 3mm Tip
- 45° Up & Down
- Bayonet



### Mucosal Dissector

- Left & Right Angle
- Spear-like Tip



### Nasal Forceps

- 13cm Working Length
- Straight & 45°Up



### Antrum Punch

- 360° Rotatable
- Back-cutting
- 10cm Working Length



### Nasal Scissors

- Pistol Grip handle
- 13cm Working Length
- Straight, Curved Right and Left



### Antrum Cannulas

- 14.5cm Total Length
- 2.5mm, 3.0mm & 4.9mm Diameter



### Kerrison Rongeurs

- 178mm & 180mm Working Length
- Angle Up & Down
- 40° Angle
- 1mm, 2mm & 3mm Biting



### One-Piece Suctions

- 5Fr, 6Fr, 7Fr, 8Fr, 9Fr, and 10Fr
- Long (150mm) & Extra Long (170mm)



Kelly Endoscopic Endonasal  
Instruments Organization Tray 1  
07-200-03MA

## Tray 2



### Decker Cup Forceps

- Oval & Round Cups
- 2mm, 2x5mm & 6mm Cup Sizes
- Straight and 45° Up
- 14.5cm & 15cmn Working Length



### Micro Scissors

- Pistol Grip Handle
- 14.5cm Working Length
- Straight and Curved Right & Left




### Spoon Forceps

- Cushing, Love-Gruenwald, Spurling Types
- 178mm Working Length
- 2x10mm, 3x10mm, 4x10mm Cup Sizes
- Straight Angle



### Dissectors

- Probe teardrop
- 2.5mm Tip
- Straight, 45°, 60° & 90°
- Straight, Up & Down Angles



### Ring Curettes

- 1.5mm, 3mm, 5mm, 7mm Cup Sizes
- 45° & 90° Up
- 45° Down



### Spatula Dissectors

- 2.4mm Tip
- Curved Up and Down



### Hemostat

- 15cm Total Length



### Mezenbaum Scissors


- 18cm Total Length
- Straight & Curved



### Mayo Scissors

- 17cm Total Length
- Straight

## Instruments Organization Trays



Kelly Endoscopic Endonasal  
Instruments Organization Tray 2  
07-200-04MA



### Weitlaner Retractor

- 13cm Total Length



### Adson Tissue Forcep

- 12cm Total Length
- 1x2mm Teeth



### Dressing Forcep

- 16cm Total Length

## Ordering Information

Product #	Description	TRAY #1	Product #	Description	TRAY #2
07-200-03MA	Kelly Endoscopic Endonasal Instruments Organization Tray 1		07-200-04MA	Kelly Endoscopic Endonasal Instruments - Tray 1 Lid & Top Tray Only	
02-034-09MA	Rhoton #3, Round Bayonetted Dissector 12cm WL, 3mm 45 Degree Up		02-202-07	Cushing, Spoon Jaw 2x10mm, 178mm Straight	
02-034-10MA	Rhoton #3, Round Bayonetted Dissector 12cm WL, 3mm 45 Degrees Down		02-202-08	Love-Gruenwald, Spoon Jaw 3x10mm, 178mm Straight	
07-181-11MA	Cottle Elevator, Double-Ended 20cm		02-202-09	Spurling, Spoon Jaw 4x10mm, 178mm Straight	
07-181-12MA	Probe, Double-Ended 19cm		07-110-10	Kelly Endonasal Yasargil Forceps, Bayonet, 23.5cm	
07-181-13MA	Nasal Forceps, 13cm WL Straight		NF-100BA	Titanium Blade Handle Rotatable Bayonet 100mm	
07-181-14MA	Nasal Forceps, 13cm WL 45 Degrees Up		02-082-09	Decker Cup Forceps, Oval 2x6mm, 15cm Straight	
07-181-15MA	Antrum Punch, 10cm WL, Backward Cutting, 360 Degrees Rotatable		02-082-20	Decker Cup Forceps, Round 2mm, 14.5cm 45 Degrees Up	
07-181-16MA	Antrum Punch, 10cm WL, Backward Cutting, 360 Degrees Rotatable		02-082-21	Decker Cup Forceps, Round 2mm, 14.5cm Straight	
07-181-17MA	Nasal Scissors, 13cm WL Straight		02-082-22	Kelly Decker Cup Forceps, Oval 2x5mm, 14.5cm 45 Degrees Up	
07-181-18MA	Kelly EndoNasal Nasal Scissors, 13cm WL Curved to Right		02-082-23	Kelly Endonasal Decker Cup Forceps, Oval 2x5mm, 14.5cm Straight	
07-181-19MA	Kelly Endonasal Nasal Scissors, 13cm WL Curved to Left		07-134-10	Kelly Endonasal Pistol Grip Micro Scissors, Straight 14.5cm WL	
07-181-20MA	Antrum Cannula w/Tear Drop Suction Control, 14.5cm Long Curved, 2.5mm Outer Dia		07-134-11	Kelly, Endonasal Pistol Grip Micro Scissors, Curved Right 14.5cm WL	
07-181-21MA	Antrum Cannula w/Tear Drop Suction Control, 14.5cm Long Curved, 3.0mm Outer Dia.		07-134-12	Kelly Endonasal Pistol Grip Micro Scissors, Curved Left 14.5cm	
07-181-22MA	Antrum Cannula w/Tear Drop Suction Control, 14.5cm Long Curved, 4.0mm Outer Dia.		06-674-00	Kelly Endoscopic Endonasal Hemostat, 15cm	
07-181-23MA	Kerrison Rongeur, 178mm 40 Degree Down, 2mm Biting		06-826-00	Mayo-Hegar Needle Holder, 16cm	
07-181-24MA	Kerrison Rongeur, 178mm 40 Degree Down, 3mm Biting		07-057-10	Weitlaner Retractor, 13cm	
07-181-25MA	Weil-Blakesley Rongeur, 13cm WL Straight 2mm		07-120-10	#3 Scalpel Handle	
07-181-26MA	Weil-Blakesley Rongeur, 13cm WL 30 Degree Up 2.6mm x 1.8mm		07-181-01	Adson Tissue Forceps, 12cm 1x2 Teeth	
07-181-27MA	Weil-Blakesley Rongeur, 2.8 x 2.5mm 13cm WL Straight		07-181-10	Dressing Forceps, 16cm	
07-181-28MA	Weil-Blakesley Rongeur, 13cm WL 30 Degree Up 3mm		09-061-10	Mayo Scissors, 17cm Straight	
02-201-01	Kerrison Rongeur, 18cm 40 Degree 1mm Up Biting		09-071-10	Metzenbaum Scissors, 18cm Straight	
02-201-02	Kerrison Rongeur, 18cm 40 Degrees 2mm Up Biting		09-071-11	Metzenbaum Scissors, 18cm Curved	
02-201-03	Kerrison Rongeur, 18cm 40 Degrees 3mm Up Biting		09-214-10	Allis Forceps, Straight 15cm	
02-034-10	Left Mucosal Dissector, Spear-Like		K-5310	Straight Blade 1mm (Box of 5)	
02-034-11	Right Mucosal Dissector, Spear-Like		K-5520	Hook Blade 30° 1mm (Box of 5)	
07-960-05L	Tapered Suction 5Fr Long		<b>Product #</b>	<b>Description</b>	<b>LIFT OUT CURETTE TRAY</b>
07-960-05XL	Tapered Suction 5Fr XL		02-034-12	Spatula Dissector, 2.4mm Curved Up	
07-960-06L	Tapered Suction 6Fr Long		02-034-13	Spatula Dissector, 2.4mm Curved Down	
07-960-06XL	Tapered Suction 6Fr XL		07-820-24	Probe Teardrop Tip, 2.5mm, Straight	
07-960-07L	Tapered Suction 7Fr Long		07-820-26	Probe Teardrop Tip, 2.5mm, 45 Degrees Up	
07-960-07XL	Tapered Suction 7Fr XL		07-820-27	Probe Teardrop Tip, 2.5mm, 90 Degrees Up	
07-960-08L	Tapered Suction 8Fr Long		07-820-29	Probe Teardrop Tip, 2.5mm, 60 Degrees Down	
07-960-08XL	Tapered Suction 8Fr XL		07-822-21	Ring Curette, 1.5mm, 45 Degrees Up	
07-960-09L	Tapered Suction 9Fr Long		07-822-23	Ring Curette, 3mm, 45 Degrees Up	
07-960-09XL	Tapered Suction 9Fr XL		07-822-25	Ring Curette, 5mm, 45 Degrees Up	
07-960-10L	Tapered Suction 10Fr Long		07-822-27	Ring Curette, 7mm, 45 Degrees Up	
07-960-10XL	Tapered Suction 10Fr XL		07-822-31	Ring Curette, 1.5mm, 90 Degrees Up	
			07-822-33	Ring Curette, 3mm, 90 Degrees Up	
			07-822-35	Ring Curette, 5mm, 90 Degrees Up	
			07-822-37	Ring Curette, 7mm, 90 Degrees Up	
			07-822-41	Ring Curette, 1.5mm, 45 Degrees Down	
			07-822-43	Ring Curette, 3mm, 45 Degrees Down	
			07-822-45	Ring Curette, 5mm, 45 Degrees Down	
			07-822-47	Ring Curette, 7mm, 45 Degrees Down	
			07-822-97	Ring Curette, 7mm, Left to Right Action	

## Notes

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---