



## Zomorodi Nano Instruments

# ZOMORODI NANO INSTRUMENTS

- Stainless steel construction for durability
- Round handles for agile manipulation and better grip
- Color-coded handles for easy identification
- Streamlined shaft for enhanced visibility


## Forceps

- 15cm, 17.5cm, 20.5cm total lengths
- 2.5cm, 5cm, or 8cm working lengths
- Carbide coated
- Straight or curved



## Scissors

- 15cm, 17.5cm, 20.5cm total lengths
- 2.5cm, 5cm, or 8cm working lengths
- Straight or curved




## Needle Holder

- 15cm, 17.5cm, 20.5cm total lengths
- 2.5cm, 5cm, or 8cm working lengths
- Carbide coated
- Straight or curved



Recommended intended use for:  
**Microsurgical procedures**



Carbide Coated Detailed Close-up

## Ordering Information

Product #	Description
07-820-00	CP and Skull Base Instrument Tray
07-828-01MA	Zomorodi Nano Forceps, Purple, Carbide Coated Tip 0.3mm, 15cm Curved, 2.5cm WL
07-828-02MA	Zomorodi Nano Forceps, Purple, Carbide Coated Tip 0.3mm, 15cm Straight, 2.5cm WL
07-828-03MA	Zomorodi Nano Scissors, Red, Sharp Tip, 15cm Curved, 2.5cm WL
07-828-04MA	Zomorodi Nano Scissors, Red, Sharp Tip, 15cm Straight, 2.5cm WL
07-828-05MA	Zomorodi Nano Forceps, Blue, Carbide Coated Tip 0.3mm, 15cm Straight, 2.5cm WL
07-828-06MA	Zomorodi Nano Forceps, Blue, Carbide Coated Tip 0.3mm, 15cm Curved, 2.5cm WL
07-828-07MA	Zomorodi Nano Forceps, Purple, Carbide Coated Tip 0.3mm, 17.5cm Curved, 5cm WL
07-828-08MA	Zomorodi Nano Forceps, Purple, Carbide Coated Tip 0.3mm, 17.5cm Straight, 5cm WL
07-828-09MA	Zomorodi Nano Scissors, Red, Curved, Sharp Tip, 5cm WL, 17.5cm TL
07-828-10MA	Zomorodi Nano Scissors, Red, Straight, Sharp Tip, 5cm WL, 17.5cm TL
07-828-11MA	Zomorodi Nano Forceps, Blue, Carbide Coated Jaw, 17.5cm Straight, 5cm WL
07-828-12MA	Zomorodi Nano Forceps, Blue, Carbide Coated Jaw, 17.5cm Straight, 5cm WL
07-828-13MA	Zomorodi Nano Forceps, Purple, Carbide Coated Tip 0.3mm, 20.5cm Curved, 8cm WL
07-828-14MA	Zomorodi Nano Forceps, Purple, Carbide Coated Tip 0.3mm, 20.5cm Straight, 8cm WL
07-828-15MA	Zomorodi Nano Scissors, Red, Curved, Sharp Tip, 8cm WL, 20.5cm TL
07-828-16MA	Zomorodi Nano Scissors, Red, Straight, Sharp Tip, 8cm WL, 20.5cm TL
07-828-17MA	Zomorodi Nano Forceps, Blue, Carbide Coated Jaw, 20.5cm Straight, 8cm WL
07-828-18MA	Zomorodi Nano Forceps, Blue, Carbide Coated Jaw, 20.5cm Curved, 8cm WL
P-071	Rubber Mat for 07-820-00 and 07-827-00 Trays

